




Cell based toxicity screening by a panel of CALUX bioassays



Dr. Peter A. Behnisch
Director
BioDetection Systems,
Amsterdam, The Netherlands
Email: peter@bds.nl



©BDS all rights reserved



Who we are



BioDetection Systems B.V. ("BDS") is a Dutch company and ISO 17025 accredited service laboratory providing biological detection systems, such as the innovative CALUX® bioassays for the determination of ultra low levels of a variety of highly potent materials.

Mission

To provide innovative bioassays and implement their use to the highest international standards.

Partner in many international projects related to water:


- EC FP6 "TECHNEAU"
- Rhine Monitoring Project (RIWA)
- EC FP7 ChemScreen for REACH/3Rs
- Dutch Projects, z. Bsp. LEOS, ZORG, Genes4Water


©BDS all rights reserved





CALUX® Principles





BDS BDS's CALUX and typical applications (1)

• DR-CALUX	(dioxins, furans, PCBs, POPs, PAKs, Ah receptor ligands)
• ERα-CALUX	(natural, synthetic and plant-like estrogens, pseudo-estrogens, DES, endocrine disrupters, tamoxifen and other anti-estrogens)
• ERβ-CALUX	(like ERα-CALUX, but more for plant-like estrogens)
• AR-CALUX	(androgens, anabolic steroids, pseudo-androgens, antibiotic growth promoters, nandrolone, flutamide e.a. steroidal drugs)
• PR-CALUX	(progesterone, MPA and other progestagene-activ compounds, steroid-based antidiuretics, aldosterone antagonists, aromatase inhibitors)

©BDS all rights reserved


BDS BDS's CALUX and typical applications (2)

• GR-CALUX	(cortisol, corticosteroids)
• TRβ-CALUX	(thyroid hormones, anti-thyroid agents, goitrogenic compounds, contrasting agents, POPs, environmental pollutants, nitrofen)
• RAR-CALUX	(retinol, vitamin A, retinoids, retinoic acid)
• NF-κB-CALUX	(endotoxins, TNFα-modulators, cytokines)
• PPARα and PPARγ-CALUX	(fatty acid derivatives, peroxisome proliferators, alkyl tin, phthalates)
• P53 CALUX	(the Guardian of DNA, genotoxic compounds)

©BDS all rights reserved



Water Hormone - Bioanalysis by CALUX



Recommendation: STOWA, the Netherlands (2010)

Endocrine disrupting chemical analyzed by ER-CALUX –
Drinking- and surface water

a) Drinking Water (RIVM) 'trigger-value': 7 ng EEQ/l.

b) Surface Water (RIVM) 'trigger-value': 1 ng EEQ/l.


In case of higher values it is recommended to find the responsible
compound and to evaluate the destruction by water treatment
plants.

Source: Mennes, W. (2004). Assessment of human health risks for estrogenic activity detected in
water samples, using the ER-CALUX assay. RIVM-notation, RIVM, Bilthoven;
see also at <http://themas.stowa.nl/Themas/Informatie.aspx?mID=7216&rd=1115&alD=1984>

©BDS all rights reserved



Dutch National Water Research (LOES 2002): Industrial and Municipal Wastewater




Temporal variation in multiple hormonal activities of surface waters located in the Dutch part of the Rhine basin





CALUX panel monitoring from Rhine surface water




Schriks et al., in prep

Appendix 2. Raw data of the present study (location Lobith), as equivalents of the given reference compound.


CALUX	ERa	PR	GR	AR	TRB
LOBITH	E2 (ng/L)	Org2058 (ng/L)	Dex (ng/L)	DHT (ng/L)	T3 (ng/L)
1-aug-07	0.029	0.031	<LOD	<LOD	<LOD
14-aug-07	0.027	0.020	<LOD	<LOD	<LOD
29-aug-07	0.026	0.027	0.92	0.017	<LOD
12-sept-07	0.73	0.039	1.7	0.034	<LOD
26-sept-07	0.04	0.038	1.4	0.025	<LOD
10-oct-07	0.032	0.028	1.4	<LOD	<LOD
24-oct-07	0.041	0.032	1.6	<LOD	<LOD
7-nov-07	0.031	0.042	2.2	0.017	<LOD
21-nov-07	0.1	0.046	2.4	0.031	<LOD
5-dec-07	0.068	0.049	2.6	0.051	<LOD
19-dec-07	0.044	0.055	1.1	<0.05	<LOD
2-jan-08	0.042	0.078	2.2	<0.05	<LOD
16-jan-08	0.075	0.068	2	<0.05	<LOD

ER-Calux: estrogen activity



- high removal rates in WTP
- care institutions: except ZE, lower activity than households
- households: RO, HA and VE high activity

GR-Calux: glucocorticoid activity



- differing removal rates in WTP
- care institutions: RO and HA high activity
- households: HA high activity



Dutch National Quality criteria's for ER CALUX (since 2002)

- Reference sample : 3 pM E2/well via Shewhart chart: 3s
- Induction of sample between LOD and EC₅₀
- Induction factor > 6
- Estradiol Standard curve linearity R² > 0.98

RIKZ-0204

Minister van Infrastructuur en Waterstaat
Rijkswaterstaat



Directoraat-Generaal
Rijksinstituut voor Kust en Zee/RIKZ

Standaardvoorschrift RIKZ/ER-CALUX® (Estrogen Receptor mediated - Chemical Activated Luciferase gene expression) bioassay.

Trefwoorden: estrogen receptor, reportergene, ER-CALUX®, estradiol

©BDS all rights reserved



General EU guideline for Screening Methods e. g. EC/1883/2006 and SANCO 10376/2011

Quantitative Analysis:

- LOD 1/5th of the relevant concentration (e.g. in Holland for surface water 1 ng EEQ/l x 1/5tel = 0,2 ng EEQ/l water)
- Standard deviations at 1/2x, 1x and 2x of the regulated levels in an acceptable range
- Participation in proficiency testing
- ISO 17025 accredited laboratory
- False Negative Rate below 5%
- Standard dilution series: linearity of standard $R^2 > 0.95$
- Triple measurement of each dilution with standard deviation below 15%; precision of 3 experiments below 30%

©BDS all rights reserved



OECD Standardization via ICCVAM/ECCVAM: REACH/ToxCast & Alternatives





.....ongoing intensive efforts of EC-ECVAM, US- ICCVAM/ToxCast and
OECD: Developments of alternative non-animals testing



©BDS all rights reserved



Evaluation for EC-CVAM: Comparison ERα CALUX® vs. other in vivo assays



©BDS all rights reserved

Sonneveld et al. (2006), Tox.Sci. 89, 173

Drinking Water Directive

⇒ Revision process is ongoing with emphasis on

- Risk assessment and management
- Microbiological parameters
- Materials in contact with drinking water
- Chemical parameters

⇒ Sound scientific knowledge, reliable and up-to-date information is of paramount importance



Tools for water quality and process control

- Electronic Nose, Electronic Tongue
- UV probe
- "Intelligent" sensors
- Bioassays





Table 1. Indicative evaluation of different effect assays for hormonal disruption (Mons, 2008).

	Sensitivity		Robustness		Time to result	Operational specifications	
	Sources	Drinking water	Operational robustness	Selectivity		Ease-of-use	Maintenance requirements
Bioassays for estrogenic activity							
ER CALUX®	4	4	5	3	3	4	3
MVLN and MELN	2	2	2	2-3	2		
T47D-KBluc	3	3			3		
YES	1	1	4	3	2-3	3	3
E-screen	4	4	4	2	1	3	3
Bioassays for androgenic activity							
AR CALUX®	3	3	3	5	3	4	3
MDA-bk2	3	3		1	3		
PALM	3	3	3	3	3		
YAS	2	1	3	3	3		
A-screen	3	3		2	1		
Bioassays for progestagenic activity							
PR CALUX®	3	3	3	5	3	4	3
TM-Luc	3	3		1	3		
Yeast based assays	2	1		4	3		
Bioassays for glucocorticoid activity							
GR CALUX®	3	3	3	5	3	4	3
TGRM-Luc	3	3		1	3		
MDA-kb2	3	3		1	3		
Bioassays for thyroid activity							
TRR CALUX®			3		3	4	3
T-Screen					1		
PC-DR-Luc							
xl58-TRE-Luc							


1 = very poor, 2 = poor, 3 = average, 4 = good, 5 = very good; when insufficient information was available, the box was left empty.

International projects and applications

- USA
- Germany
- Japan
- European Union




Estrogen-like compounds by ER CALUX





Androgen-like compounds by AR CALUX





Take home message

- Multiple biotectors or Effect based analysis tools have been evaluated in many international projects for many environmental applications – they are currently also in several countries parallel in the last evaluation phase (EU, China, Japan, Australia)
- Endocrine disrupters are not only female hormones – male and other hormones (PR, GR, PPAR, RXR, TR) need more attention!
- Complex mixture cocktails and multi-pollutants effects are relevant and can be now evaluated by a panel of CALUX tests
- Search for R&D partners and policymakers to move Multiple biotectors or Effect based analysis tools forward..interested ?

©BDS all rights reserved



**6th BioDetectors Workshop
in Amsterdam in May 2012**

