

**Quality Assurance Project Plan (QAPP) for the State Water Resources Control  
Board *Ceriodaphnia dubia* Study - Baseline Interlaboratory Testing**

**08/08/22**

## Title and Approval Sheet

Title: Quality Assurance Project Plan (QAPP) for the State Water Resources Control Board *Ceriodaphnia dubia* Interlaboratory Comparison Study

Prepared by: \_\_\_\_\_

QAPP Approval Date: \_\_\_\_\_

---

Approved:	Organization & Role	Name	Date
Approved:	Organization & Role	Name	Date
Approved:	Organization & Role	Name	Date
Approved:	Organization & Role	Name	Date
Approved:	Organization & Role	Name	Date
Approved:	Organization & Role	Name	Date

# Table of Contents

SECTION	PAGE
<b>1 PROJECT MANAGEMENT</b> .....	<b>1</b>
1.1 Title and Approval Sheet (QA/R-5 Element A1).....	1
1.2 Table of Contents (QA/R-5 Element A2).....	1
1.3 Distribution List (QA/R-5 Element A3).....	1
1.4 Introduction and Project Organization (QA/R-5 Element A4).....	1
1.5 Problem Definition and Background (QA/R-5 Element A5) .....	2
1.6 Project Description (QA/R-5 Element A6).....	5
1.7 Quality Objectives and Criteria (QA/R-5 Element A7).....	10
1.8 Training and Certification Requirements (QA/R-5 Element A8).....	16
1.9 Documentation and Records (QA/R-5 Element A9).....	16
<b>2 DATA GENERATION AND ACQUISITION</b> .....	<b>16</b>
2.1 Experimental Design (QA/R-5 Element B1) .....	16
2.2 Sampling Methods (QA/R-5 Element B2) .....	17
2.3 Sample Handling and Custody (QA/R-5 Element B3).....	18
2.4 Analytical Methods (QA/R-5 Element B4) .....	18
2.5 Quality Control (QA/R-5 Element B5) .....	19
2.6 Instrument/Equipment Testing, Inspection, and Maintenance (QA/R-5 Element B6).....	21
2.7 Instrument/Equipment Calibration and Frequency (QA/R-5 Element B7).....	21
2.8 Inspection/Acceptance of Supplies and Consumables (QA/R-5 Element B8).....	22
2.9 Non-direct Measurements (QA/R-5 Element B9).....	23
2.10 Data Management (QA/R-5 Element B10).....	23
<b>3 ASSESSMENT AND OVERSIGHT</b> .....	<b>23</b>
3.1 Assessment and Response Actions (QA/R-5 Element C1) .....	23
3.2 Reporting (QA/R-5 Element C2) .....	23
<b>4 DATA REVIEW AND USABILITY</b> .....	<b>24</b>
4.1 Data Review, Verification, and Validation (QA/R-5 Element D1).....	24
4.2 Verification and Validation Methods (QA/R-5 Element D3) .....	24
4.3 Reconciliation with User Requirements (QA/R-5 Element D3) .....	26
<b>5 REFERENCES</b> .....	<b>27</b>

## LIST OF TABLES

- Table 1. Interlaboratory study baseline tests to be conducted by participating laboratories.
- Table 2. Summary of test conditions and acceptability criteria for conducting chronic *C. dubia* toxicity tests.
- Table 3. Water quality measurements and methods.
- Table 4. Chemical measurements and methods for aqueous samples.
- Table 5. Laboratory quality control samples for the *C. dubia* interlaboratory toxicity tests
- Table 6. Quality control samples for chemical analyses.

## APPENDICES

- Appendix A - Laboratory Contact Information
- Appendix B – Data Quality Objectives
- Appendix C – Toxicity Data Submittal Instructions
- Appendix D – Study Plan and Logistics

## Distribution List

<b>SCCWRP Co-Project Manager/ Technical Coordinator</b>	Alvina Mehinto
<b>SCCWRP Co-Project Manager</b>	Ken Schiff
<b>SCCWRP Quality Assurance Officer</b>	Darrin Greenstein
<b>SCCWRP Biostatistician</b>	David Gillett and Jing Zhang
<b>SCCWRP Data Manager</b>	Paul Smith
<b>Laboratory 1 Project Manager</b>	_____
	Name
<b>Laboratory 1 QA Manager</b>	_____
	Name
<b>Laboratory 2 Project Manager</b>	_____
	Name
<b>Laboratory 2 QA Manager</b>	_____
	Name
<b>Laboratory 3 Project Manager</b>	_____
	Name
<b>Laboratory 3 QA Manager</b>	_____
	Name
<b>Laboratory 4 Project Manager</b>	_____
	Name
<b>Laboratory 4 QA Manager</b>	_____
	Name
<b>Laboratory 5 Project Manager</b>	_____
	Name
<b>Laboratory 5 QA Manager</b>	_____
	Name
<b>Laboratory 6 Project Manager</b>	_____
	Name
<b>Laboratory 6 QA Manager</b>	_____
	Name

**Laboratory 7 Project Manager**

Name

**Laboratory 7 QA Manager**

Name

**Laboratory 8 Project Manager**

Name

**Laboratory 8 QA Manager**

Name

**Laboratory 9 Project Manager**

Name

**Laboratory 9 QA Manager**

Name

**Laboratory 10 Project Manager**

Name

**Laboratory 10 QA Manager**

Name

**Laboratory 11 Project Manager**

Name

**Laboratory 11 QA Manager**

Name

**Laboratory 12 Project Manager**

Name

**Laboratory 12 QA Manager**

Name

# 1 PROJECT MANAGEMENT

## 1.1 TITLE AND APPROVAL SHEET (QA/R-5 ELEMENT A1)

See page i.

## 1.2 TABLE OF CONTENTS (QA/R-5 ELEMENT A2)

See page ii.

## 1.3 DISTRIBUTION LIST (QA/R-5 ELEMENT A3)

See page iv.

## 1.4 INTRODUCTION AND PROJECT ORGANIZATION (QA/R-5 ELEMENT A4)

This Quality Assurance Project Plan (QAPP) was developed according to United States Environmental Protection Agency (EPA 2002a, 2016) guidance to inform the implementation of the Interlaboratory Study (ILS) component of the California State Water Resources Control Board (SWB) study titled *Development of Quality Assurance Recommendations for the Ceriodaphnia Toxicity Test* (SWB *C. dubia* study). The ILS was informed by SCCWRP discussions with an Expert Science Panel (ESP) and a Stakeholder Advisory Committee (SAC), ASTM (1996) *Standard Practice for Conducting an Interlaboratory Study to Determine the Precision of a Test Method*, and the EPA (2001a,b) *Interlaboratory Variability Study of EPA Short-term Chronic and Acute Whole Effluent Toxicity Test Methods*.

Roles and responsibilities for individuals needed to implement task planning, management and oversight, sampling/sample preparation, sample analysis, and data management for the ILS are summarized below.

The Southern California Coastal Watershed Research Project (SCCWRP) is implementing this work on behalf of the SWB in coordination with commercial and municipal laboratories.

**Co-Project Manager /Technical Coordinator** – Dr. Alvina Mehinto is the SCCWRP project manager responsible for overall project approval and the lead coordinator responsible for the distribution of this QAPP, providing updates to all participants and corresponding with the Stakeholder Advisory Committee and Expert Science Panel.

**Co-Project Manager** – Ken Schiff is responsible for SCCWRP’s oversight of the project to ensure that the needed resources are available to complete this project.

**Quality Assurance (QA) Officer** – Darrin Greenstein will oversee the preparation and shipping of all samples and conduct quality assurance oversight of the tests conducted by laboratories and chemical analyses performed by the analytical laboratory to ensure submitted data are accurately reported and meet quality criteria described in this QAPP.

**Biostatistician** – Drs David Gillett (SCCWRP) and Jing Zhang (Miami University) are the biostatisticians knowledgeable regarding statistical procedures needed to meet the goals of this study and with the materials being tested in order to ensure that data to be collected will meet the statistical needs to achieve the study goals. The biostatisticians will also assist in interpreting the results of data analysis.

**Data Manager** – Paul Smith will coordinate data compilation and will perform calculations to evaluate data.

Sixteen laboratories participated in the SCCWRP's historical data and SOPs request. Among them, twelve (12) laboratories have committed to participating in the ILS (**Appendix A**). These laboratories represent most Environmental Laboratory Accreditation Program (ELAP) accredited laboratories; thus, they are a reasonable representation of qualified laboratories applicable to toxicity testing in California (i.e., they are ELAP accredited for the test method) and this subset can adequately characterize test reproducibility even if a few do not participate. The ILS is not restricted to a group of laboratories judged to be exceptionally qualified and equipped for the ILS (ASTM 1996).

**Laboratory Coordinators** – Each participating laboratory will identify a laboratory coordinator to oversee their implementation of testing associated with this ILS. They will also coordinate communication between the laboratory and the SCCWRP Technical Coordinator and Quality Assurance Officer.

**Laboratory QA Manager** – Each participating laboratory will identify a QA manager to ensure that testing is conducted in accordance with laboratory SOPs and testing protocols described in this QAPP. They also perform QA reviews of submitted data to confirm that results are accurately reported and meet quality criteria described in this QAPP.

## 1.5 PROBLEM DEFINITION AND BACKGROUND (QA/R-5 ELEMENT A5)

This section describes the rationale for this project and relevant background information.



### 1.5.1 SWB *C. dubia* Study

The SWB adopted Resolution No. 2020-0044<sup>1</sup> in 2020 to establish the *Water Quality Control Plan for Inland Surface Waters, Enclosed Bays, and Estuaries of California and Adopting Toxicity Provisions* (Toxicity Provisions). Public comments at the SWB adoption hearing on December 1, 2020, communicated that dischargers were primarily concerned with interlaboratory variability for the chronic *C. dubia* method (achieving the same test result among multiple laboratories testing the same sample) and incorrect indications of toxicity on non-toxic samples.

To address stakeholders' concerns, the SWB is funding a *C. dubia* QA evaluation study, facilitated by SCCWRP, that aims to provide laboratory technique guidance to: (a) improve the consistency of the execution of the *C. dubia* test method within each testing laboratory; and (b) improve the consistency and comparability of *C. dubia* test results among testing laboratories, while retaining the necessary flexibility for environmental relevance. The ILS represents one of the tasks being conducted as part of the SWB *C. dubia* study.

The SWB *C. dubia* study is divided into four primary tasks. These are described in the Conceptual Workplan (SCCWRP 2021b) and summarized below. The baseline ILS described by this QAPP is Task 2.3. The QAPP describes the approach, overall methodology, and logistics that will be used to conduct the study

- Task 1 – Establish a governance structure for the study (completed).
- Task 2.1 and 2.2 – Inventory of protocols and historical control data for the *C. dubia* chronic WET test from all ELAP accredited laboratories (completed).
- Task 2.3 (optional) – A baseline ILS will be conducted using split samples to assess inter- and intra-laboratory variability (described in this QAPP). The split sample analysis data will supplement historical data analyses and confirm possible sources of test variation.
- Task 3 – Targeted experiments will be conducted among a few laboratories to standardize select test parameters to minimize inter- and intra-laboratory variability in test results.
- Task 4 – A second round of split-sample interlaboratory testing will be conducted to determine if recommendations developed in Task 3 are successful in reducing inter- and/or intra-laboratory variability among a wide range of laboratories.
- Task 5 – A final report will be prepared that describes study results and recommendations for reducing inter- and intra-laboratory variability.

---

<sup>1</sup> In 2021, the State Water Board adopted Resolution No. 2021-0044 to rescind the December 1, 2020, establishment of water quality control plan for inland surface waters, enclosed bays, and estuaries of California and confirmed that the “Toxicity Provisions” were adopted as state policy for water quality control for all waters of the state.

### 1.5.2 *C. dubia* Test Variability

The *C. dubia* short-term chronic survival and reproduction toxicity test (EPA 2002b, 2016, test method 1002.0) is an established whole effluent toxicity (WET) test commonly used in regulatory and monitoring programs, including the Toxicity Provisions (2020b). EPA conducted two interlaboratory studies associated with the development of the *C. dubia* chronic toxicity test. In the first national study of interlaboratory precision performed in 1991, the coefficient of variation (CV), a statistical measure of the relative dispersion of data points in a data series around the mean, was 72.9% for the 25 percent inhibition concentration (IC25) determined by 155 participating laboratories (EPA 2002b). Over time, this precision improved and in 2000 the IC25 from 35 interlaboratory chronic toxicity tests with *C. dubia* had a CV of 35% in the combined dataset consisting of valid tests with KCl spiked effluent and KCl spiked receiving water (EPA 2002b, 2016, 2001a).

EPA (2000) describes toxicity test variability and discusses factors that can contribute to inter- and intra-laboratory variability as follows.

*Variability is inherent in any analytical procedure. The precision of a method describes the closeness of agreement between test results obtained from repeated testing of a prescribed method. WET test precision can be categorized by: 1) intratest (within-test) variability, 2) intralaboratory (within-laboratory) variability, and 3) interlaboratory (between-laboratory) variability. Intratest variability can be attributed to variables such as the number of treatment replicates, the number of test organisms exposed per replicate, and the sensitivity differences between individual organisms (i.e., genetic variability). Intralaboratory variability is that which is measured when tests are conducted under reasonably constant conditions in the same laboratory (e.g., reference toxicant or effluent sample tested over time). Sources of intralaboratory variability include those factors described for intratest variability, as well as differences: 1) in test conditions (e.g., seasonal differences in dilution water quality, differences in environmental conditions), 2) from test to test in organism condition/health, and 3) in analyst performance from test to test. Interlaboratory variability reflects the degree of precision that is measured when the same sample or reference toxicant is analyzed by multiple laboratories using the same methods. Variability measured between laboratories is a consequence of variability associated with both intra-test and intra-laboratory variability factors, as well as differences allowed within the test methods themselves (e.g., source of dilution water), technician training programs, sample and organism culturing/shipping effects, testing protocols, food quality, and testing facilities.*

### 1.5.3 Historical Data Evaluation

SCCWRP (2022) conducted a historical data analysis (SWB *C. dubia* study Task 2.2) intended to "... use test results to identify a subset of factors for which there is a certain allowable amount of flexibility in the U.S. EPA standard method (e.g., age at start, test duration, ionic composition, water chemistry, culture feeding regime, and others) that appear to influence the results between

laboratories.” This evaluation of control and reference toxicant data collected over a 2 to 5-year period, from 16 laboratories accredited by ELAP to perform the chronic *C. dubia* toxicity test did not identify any single consistent factor associated with inter- or intra-laboratory variability. As the result the ESP could not identify a subset of laboratory techniques that should be further optimized. Instead, the ESP and SAC recommended proceeding with an interlaboratory split-sample comparison exercise (SWB *C. dubia* study Task 2.3) to collect additional data to support the SWB *C. dubia* study given the lack of any clear driver of variability in the historical data.

## **1.6 PROJECT DESCRIPTION (QA/R-5 ELEMENT A6)**

This section summarizes the work to be performed and any decision(s) to be made or outcomes expected from information to be obtained.

### **1.6.1 Study Goals and Objectives**

The objective of this study is to build on previous efforts and investigate a variety of possible sources of variability in the *C. dubia* reproduction test conducted by ELAP-accredited laboratories. The primary goals of the SWB *C. dubia* study, in collaboration with stakeholders and laboratories, are as follows.

- 1) Evaluate laboratory performance among those accredited by the state of California (i.e., ELAP).
- 2) Investigate factors that can lead to inter- and intra-laboratory test variability and decrease confidence in assessments of toxicity.
- 3) Recommend specific guidance for laboratory techniques to improve laboratory performance reduce intra- and inter- laboratory variability (SWB 2020a, SCCWRP 2021a, 2022).

The ILS is intended to address the following two general questions from the ESP and SAC.

- 1) Which lab practice(s) should be standardized to reduce inter- and intra-laboratory variability?
- 2) Does standardizing lab practices improve consistency and comparability in *C. dubia* test results?

### **1.6.2 Approach**

The baseline ILS testing is summarized in this section to provide an overview and rationale for the planned approach.

The study will assess the toxicity of different samples provided by SCCWRP using the chronic *C. dubia* WET test (EPA 2002b). The study design, developed by the ESP, focuses on *C. dubia* reproductive endpoints (i.e., 25 and 50% inhibition concentrations, IC 25 and IC50). Test samples were selected to assess interlaboratory agreement in estimated IC values and the impact of different dilution water recipes on test outcome.

Twelve out of the 16 laboratories that provided historical data, are participating in the study. Among the laboratories missing, two are no longer ELAP accredited, two have scheduling conflicts and one laboratory did not respond to our request.

### *Samples*

The three sets of samples to be tested as part of this ILS are described below and in

#### **Table 1.**

- Sample 1: Moderately hard dilution water recipe #1 (EPA MH) to be tested at full strength (i.e., 100%). This sample shall be tested along with one (1) laboratory control consisting of the lab's own dilution water recipe.

Sample 1 will be prepared using EPA MH water following the protocol described in the EPA manual. This recipe is only used by two (2) of the 12 participating labs. The other laboratories are either supplementing the EPA MH water with selenium and/or vitamins in various quantities, using the modified EPA method or using mineral water. Due to the issues documented with the use of EPA MH alone by several laboratories (i.e., low reproduction and poor culture viability), this recipe was not selected for the dilution series. This sample will be used to compare intra- and inter-laboratory performance using the EPA recommended method.

- Sample 2A: Moderately hard dilution water recipe #2 (Perrier) to be tested at full strength (i.e., 100%). This sample shall be tested along with one (1) laboratory control consisting of the lab's own dilution water recipe.
- Sample 2B-F: 5 concentrations of sodium chloride (NaCl; 2000, 1000, 500, 250, 125 mg/L) diluted in moderately hard water recipe #2, Perrier water. All dilutions will be prepared at SCCWRP. Therefore, the samples shall be tested as is (i.e., no additional sample dilution allowed) along with one (1) laboratory control consisting of the lab's own dilution water recipe.

Sample 2A-F will use Perrier water, a recipe used by four (4) of 12 participating laboratories to prepare the dilution series with NaCl. The concentrations of NaCl were determined based on previous work from the EPA and historical data from the participating laboratories using this salt as their reference toxicant. Sample 2 will be run

with two laboratory controls (i.e., sample 2A + one laboratory control; and sample 2B-F + one laboratory control) to ensure that individual laboratories generate sufficient control data to assess intra-laboratory variability. However, laboratories will treat sample 2A-F as one test to be conducted simultaneously using blocked randomization from the same known parentage. Sample 2 will be used to assess interlaboratory agreement for *C. dubia* reproductive endpoints IC25 and IC50.

- Sample 3: NaCl will be provided (7 g as a solid) to each lab with instructions to prepare 5 dilutions using the lab’s own dilution water. This serial dilution will be tested along with one (1) laboratory control consisting of the lab’s own dilution water recipe. Note that

Participating laboratories are currently using different types and concentrations of reference toxicant, limiting our ability for inter-laboratory comparison across a large number of laboratories. Here Sample 3 will replace the reference toxicant and laboratories will prepare a NaCl (i.e., 0, 125, 250, 500, 1000 and 2000 mg/L) using their own control water. Data generated will be used to evaluate the impact of dilution water and lab techniques on reproductive endpoints.

**Table 1. Interlaboratory study baseline tests to be conducted by participating laboratories.**

Sample Type	Number of Samples	Dilution Series	Number of Rounds	Number of Labs
Sample 1- dilution water recipe #1 (EPA MH)	1	No	3	12
Sample 2A- dilution water recipe #2 (Perrier)	1	No	3	12
Sample 2B-F- NaCl in dilution water recipe #2 (Perrier)	5	No	3	12
Sample 3- NaCl solid sample (to be diluted by the labs using their own dilution water)	1	Yes	3	12

Baseline ILS testing will be conducted with three rounds of tests (i.e., repeat testing with newly prepared samples in each of three rounds of testing). This approach will yield up to 144 test results (4 tests per batch x 3 batches x 12 participating laboratories), as shown in **Table 1** above.

### *Sample Preparation and Testing*

Sample preparation and testing will be conducted as follows.

- New samples will be prepared by SCCWRP for each round of testing and allowed to equilibrate for 48 hours.
- Prior to starting the first round of testing, SCCWRP will prepare bulk samples that will be tested by one laboratory to verify that their level of toxicity is within the expected range. If unspiked samples are not toxic and a dose-response is observed, a new batch of

samples will then be prepared, split in cubitainers, and shipped to the laboratories. Since all methods and equipment will be the same for subsequent rounds, this preliminary testing will only be carried out for round one.

- The water samples will be split into cubitainers for distribution to the laboratories. The original bulk sample and each split sample will be analyzed for alkalinity, conductivity, hardness, and pH prior to shipment to document consistency with the parent sample and among the splits. A subset of the split-samples will also be collected by SCCWRP for ion composition analyses.
- Samples will be maintained at  $4^{\circ}\text{C} \pm 2$  prior to use.
- Laboratory testing will be conducted promptly (i.e., within 24 hours of receiving the split samples) so that sample holding times are not exceeded.
- The ILS will be conducted over an eight-week period.
- Chronic *C. dubia* testing will include:
  - Static exposures with daily renewal using the single cubitainer supplied by SCCWRP.
  - 10 replicate chambers per sample/dilution concentration.
  - Survival and reproduction (i.e., neonates per female) will be documented daily.
  - Each sample will be tested with a separate laboratory control (i.e., the standard control water used by each laboratory). One control will be used for samples tested as a dilution series.
  - Water chemistry measurements will include ionic composition (see section 1.7.3). SCCWRP will collect samples of the test solutions before they are shipped to the laboratories. Laboratories will collect a sample of their standard control water at test initiation using the containers provided by SCCWRP and ship them to SCCWRP. All samples for ionic composition will be measured by Physis.
  - Test set-up must use the randomized blocking by known parentage, using only brood board chambers with a minimum of 8 neonates from the adult on test initiation day. Each test will be treated as independent for blocking randomization except for sample 2A and 2B-F and the associated laboratory controls that must be blocked from the same known parentage.
  - Tests will be conducted over 8-days (i.e., carried out to 192 hours).

### *Laboratory Documentation*

The following information will be collected and reported by each laboratory in electronic format.

- Alkalinity, conductivity, dissolved oxygen (DO), hardness, pH, and temperature will be measured and recorded for each sample upon receipt by each laboratory.

- Counts of males, unhealthy and dead adults, and dead neonates will be documented in the brood boards used to set up the tests.
- Specific beginning and end time window for age of neonates at test initiation will be recorded.
- Test solutions will be renewed daily at 24 hours, within a +/- 1 hour window, to enhance comparability of neonate counts among laboratories. The specific time of renewal (hours and minutes) shall be recorded.
- Water quality parameters (dissolved oxygen, air and water temperature, pH, and conductivity,) will be measured before and after daily renewal, following the laboratories' routine procedures. The laboratories will document their procedures including methods used for sample collection, analysis (including measurement devices), and volume of sample needed. It is also recommended that the laboratories measure water temperature continuously.
- Dilution water recipe, the food recipe, and preparation dates will be reported.
- Light intensity will be measured once at test initiation within the location of each test.
- Air temperature will be measured twice daily at the location of the testing chambers.
- Control charts for reference toxicant testing over the last 12 months (i.e., October 2021 to September 2022) will be provided to SCCWRP. For laboratories conducting multiple reference toxicant tests per month, the first test of the month will be provided. Charts will include IC50 and IC25 data and a tabular listing of the individual IC values.
- Laboratories will use the bench sheets templates provided by SCCWRP and submit scanned copies to SCCWRP as part of their data submission package. The laboratories may use their own bench sheets as long as they include all data categories found in the SCCWRP templates.

### **1.6.3 Possible Outcomes and Decisions**

The ILS will produce a baseline dataset that can be evaluated by SCCWRP and discussed by the SAC and ESP to focus subsequent tasks in the SWB *C. dubia* study. A possible outcome is that baseline ILS data can be used to identify which factors or lab techniques are associated with intra- and/or inter-laboratory variability and can be further investigated in Tasks 3 and 4 of the SWB *C. dubia* study. These investigations will lead to recommendations to the SWB to reduce sources of variability and improve consistency and comparability in *C. dubia* chronic toxicity test results. The SWB will make decisions regarding how to use any recommendations developed from this study.

## 1.7 QUALITY OBJECTIVES AND CRITERIA (QA/R-5 ELEMENT A7)

This section describes how data of known and acceptable quality will be developed to satisfy the study objectives.

### 1.7.1 Data Quality Objectives

A systematic project planning process, the seven-step Data Quality Objectives (DQOs) process, was used to determine study boundaries and specifications needed to develop a study that will support the qualitative and quantitative derivation of data, based on EPA (2006) guidance. DQOs are based on a seven-step process.

- Step 1: State the Problem
- Step 2: Identify the Goal of the Study
- Step 3: Identify Information Inputs
- Step 4: Define the Boundaries of the Study
- Step 5: Develop the Analytical Approach
- Step 6: Specify Performance or Acceptance Criteria
- Step 7: Develop the Plan for Obtaining Data

DQOs specify the data types, quality, quantity, and uses, as needed to make decisions, and are the basis for designing data and sample collection activities. These minimum standards are developed to provide data of known quality to address the study objectives. **Appendix B** presents DQOs for this ILS.

### 1.7.2 Data Quality Indicators

Performance and acceptance criteria are described in terms of Data Quality Indicators (DQIs). DQIs for this project consist of qualitative and quantitative indicators, including precision, accuracy, representativeness, completeness, and comparability.

#### *Precision*

Tests performed on presumably identical materials in presumably identical circumstances do not typically yield identical results. An indication of a test method's consistency is its precision – a measure of agreement for repeated analyses under similar conditions over a relatively short period of time. The CV is a simple statistic to describe precision by comparing the results of replicate analyses within a single lab (i.e., intra-laboratory precision or *repeatability*) or among different laboratories (i.e., inter-laboratory precision or *reproducibility*).



Coefficient of Variation (CV) = standard deviation / mean

The use of a CV removes units from the measurement and allows comparisons of variability among different types of test methods or toxicity endpoints (i.e., neonates per female and analytical chemistry). However, there is no common agreement on an acceptable standard deviation (EPA 2000).

The relative percent difference (RPD) is another measure of inter- or intra-laboratory variability among repeated measurements.

Relative Percent Difference (RPD) =  $((M - MD) / \text{mean}) * 100$

Where:

M: Measured value associated with the primary sample

MD: Measured value associated with the duplicate or replicate analysis

Mean: the mean of the two measurements  $((M + MD)/2)$ .

Precision estimates should be obtained through the efforts of qualified laboratories and personnel operating under prevailing conditions when the test method is used in practice (ASTM 1996).

### *Accuracy*

Accuracy is the degree of agreement between a measured value and the “true” or expected value. As such, it represents an estimate of total error from a single measurement, including both systematic error (“bias”) and random error that may reflect variability because of imprecision. Accuracy is typically expressed in terms of percent recoveries determined from results of matrix spike/matrix spike duplicates, laboratory control standards, and ongoing evaluation of calibration verification information. Such measurements are appropriate for chemical analyses but not for toxicity testing, which is method dependent and, thus, is subject to the variability inherent among test organisms.

### *Representativeness*

Representativeness is the degree to which sample data accurately expresses the characteristics of a population of samples, parameter variations at a sampling point, or an environmental condition. It is a qualitative parameter achieved through proper sampling program design and use of appropriate sampling strategies and techniques. For example, holding time requirements are intended to ensure the representativeness of conditions at the time the samples are collected. The use of quality control (QC) samples that are similar in composition to samples being measured provides a means of estimating precision and accuracy that are representative of sample

measurements. In addition to sample holding times, a factor that affects representativeness in the ILS study is split-sample homogeneity. This will be evaluated by comparing water quality measurements among split samples and comparisons with the parent sample.

### *Completeness*

Completeness can be defined both qualitatively and quantitatively. Qualitative completeness is determined as a function of all factors that contribute to sampling. Quantitative completeness is calculated as the percentage of measurements that are judged to be valid compared to the total number of measurements planned. Effectively, it measures the amount of data available for valid measurements compared to the amount lost or destroyed, and it is strictly defined as the ratio of the number of usable data points over the possible number of data points, by method/matrix. A completeness goal for the baseline ILS is 100 percent.

### *Comparability*

Comparability is a qualitative indicator of the confidence with which one data set can be compared to another. Confidence is achieved by maintaining standard techniques and procedures for collecting and analyzing representative samples and reporting the analytical results in standard units. ELAP certified laboratory analytical and toxicity testing and EPA approved methods will be used to ensure comparable data are produced.

### **1.7.3 Measurement Quality Objectives**

Laboratory QC results must meet the performance/acceptance limits detailed in the applicable measurement quality objectives (MQOs), or acceptance thresholds and goals, for chemical analyses and toxicity testing. These specify how much error is acceptable and are often expressed in terms of DQIs. Chronic *C. dubia* toxicity test conditions and test acceptability criteria (TAC) from EPA (2002b) are summarized in **Table 2**. DQIs for water quality measurements and methods are presented in **Table 3**.

MQOs for chemical analyses conducted by the analytical laboratory, Physis, are presented in

### **Table 4.**

**Table 2. Summary of test conditions and acceptability criteria for conducting chronic *C. dubia* toxicity tests.**

<b>Parameter</b>	<b>Measurement Quality Objective</b>
Test type	Static renewal (required)
Test temperature (°C)	25.0 ± 1.0°C (recommended) Test temperatures must not deviate (i.e., maximum minus minimum temperature) by more than 3.0°C during the test (required)
Light quality	Ambient laboratory illumination (recommended)
Light intensity	10-20 µE/m <sup>2</sup> /s or 50-100 ft-c (ambient laboratory levels) (recommended)
Photoperiod	16 h light, 8 h dark (recommended)
Test chamber size	30 mL (recommended minimum)
Test solution volume	15 mL (recommended minimum)
Renewals	Daily (required)
Age of test organisms	Less than 24 h; and all released within an 8-h period (required)
Number of organisms per replicate	One (1) per replicate. Assigned using blocking by known parentage (required)
Number of replicate test chambers per concentration	10 (required)
Feeding regime	Feed 0.1 mL each of YCT and algal suspension per test chamber daily (recommended)
Cleaning	Use freshly cleaned glass beakers or new plastic cups daily (recommended)
Control water	Uncontaminated source of receiving or other natural water, synthetic water prepared using MILLIPORE MILLI-Q® or equivalent deionized water and reagent grade chemicals or DMW
Test Concentrations	All split water samples provided for this ILS will be tested without dilution by the laboratories.
Dilution Factor	Sample 3 (provided as solid) will be tested in 5 dilutions (2000, 1000, 500, 250 and 125 mg/L)
Test Duration	8 days (required)* Note that EPA (2002a) defines the test duration as “Until 60% or more of surviving control females have three broods (maximum test duration 8 days) (required)”
Endpoints	Survival and reproduction (required)
Test acceptability criteria	80% or greater survival of all control organisms and an average of 15 or more live young per surviving female in the control solutions. 60% of surviving control females must produce three broods (required)
Sampling requirements	A single split sample will be provided for testing.
Sample volume required	One (1) L/day (recommended)
Water quality measurements	Initial: DO, pH, temperature, conductivity, alkalinity, and hardness (required)

**Table 2. Summary of test conditions and acceptability criteria for conducting chronic *C. dubia* toxicity tests.**

Parameter	Measurement Quality Objective
	Daily before and after renewal: DO, pH, air and water temperature, and conductivity (required); continuous water temperature measure (recommended)
	Final: DO, pH, temperature, conductivity, alkalinity, and hardness (required)
Aeration	None; unless DO falls below 4.0 mg/L or the sample has high likelihood to have substantial BOD
<p>Notes:</p> <p>Source: EPA 2002b (Method 1002.0)</p> <p>Each test condition listed above is identified as required or recommended for the purposes of reviewing test data submitted under NPDES permits. These requirements are applicable to the routine test procedures that would apply under the Toxicity Provisions and will be followed for the ILS.</p> <p>*ILS specific requirement</p> <p>BOD – Biological Oxygen Demand</p> <p>DO – dissolved oxygen</p> <p>DMW – demineralized water</p> <p>ILS – Interlaboratory study</p> <p>NPDES – National Pollutant Discharge Elimination System</p> <p>EPA – United States Environmental Protection Agency</p> <p>YCT – Yeast, Cerophyl®, trout chow</p>	

**Table 3. Water quality measurements and methods to be used by the participating testing laboratories.**

Parameter	Method <sup>1</sup>	Units	Accuracy	Precision <sup>2</sup>	Resolution <sup>3</sup>
Dissolved Oxygen	Meter or probe	mg/L	± 0.5	± 0.5 or 10%	± 0.1
pH	Meter or probe	pH units	± 0.2	± 0.2	± 0.1
Conductivity	Meter or probe	µS/cm @25°C	± 2	± 2 or 10%	± 1
Temperature	Meter, probe, or HOBO	°C	± 1	± 1 or 10%	± 0.1
Alkalinity	Titration or spectrophotometric	mg/L CaCO3	± 2	± 2 or 10%	± 1
Hardness	Titration or spectrophotometric	mg/L CaCO3	± 2	± 2 or 10%	± 1

Notes:

Source: Conventional Parameters in Fresh and Marine Water

([https://www.waterboards.ca.gov/water\\_issues/programs/swamp/docs/mqo/conventional\\_parameters\\_in\\_fresh\\_water\\_and\\_marine\\_water.pdf](https://www.waterboards.ca.gov/water_issues/programs/swamp/docs/mqo/conventional_parameters_in_fresh_water_and_marine_water.pdf))

Notes: Readings should be recorded with significant figures as shown in the resolution column.

<sup>1</sup> Method used must be documented in the data submission form.

<sup>2</sup> Precision is defined by the Relative Percent Difference (RPD). This is the difference between two repeated measurements expressed as a percentage. %RPD = (sample result – duplicate result) / sample result x 100.

<sup>3</sup> Resolution refers to the capability of a method or instrument to recognize small differences between values. This term is often used to assess if an instrument or method is useful to a study and is provided by the manufacturer.

**Table 4. Chemical measurements and methods for analyses of aqueous samples by the analytical lab.<sup>1</sup>**

Parameter	Method	Units	Method Detection Limit	Reporting Limit	Holding time <sup>2</sup>
Alkalinity, Total	SM2320B	mg/L	1.0	5.0	
Ammonia, Total	SM 4500-NH3	mg/L	0.03	1.0	
Hardness, Total	EPA 200.7	mg/L	0.19	1.0	
Ions					
Carbonate	SM 2320 B	mg/L	1	5	14 days
Bicarbonate	SM 2320 B	mg/L	1	5	14 days
Chloride	EPA 300.0	mg/L	0.01	0.05	28 days
.....Fluoride	EPA 300.0	mg/L	0.01	0.05	28 days
.....Nitrate	EPA 300.0	mg/L	0.01	0.05	28 days
Sulfate	EPA 300.0	mg/L	0.01	0.05	28 days
Major Cations					
Calcium	EPA 200.8	mg/L	0.05	0.1	180 days
Magnesium	EPA 200.8	mg/L	0.05	0.1	180 days
Potassium	EPA 200.8	mg/L	5	10	180 days
Sodium	EPA 200.8	mg/L	5	10	180 days
Selenium, Total	EPA 200.8	µg/L	0.02	0.07	180 days

<sup>1</sup> All 500-mL samples collected by the laboratories (i.e., their dilution water) and SCCWRP (i.e., test samples) will be analyzed for the suite of parameters listed in this table.

<sup>2</sup>All samples will be held at 4°C in the dark until analyzed. No additional preservation is needed.

## 1.8 TRAINING AND CERTIFICATION REQUIREMENTS (QA/R-5 ELEMENT A8)

Baseline ILS participants will have appropriate training and experience conducting their tasks, including conducting the short-term chronic *C. dubia* toxicity test, data collection, data analysis, and reporting. Toxicity testing laboratory staff will be current on all training protocols required for implementing the *C. dubia* test, conducting water quality measurements, and record keeping. These activities will be performed according to their own laboratory SOPs, unless otherwise required to meet project specific requirements, such as additional data collection required by this QAPP (e.g., extending the test duration through day-8 may differ from standard laboratory practices, counting live and dead young, etc.). Analytical and toxicity testing laboratories will be currently or recently accredited by ELAP to conduct all testing procedures they are assigned to perform.

## 1.9 DOCUMENTATION AND RECORDS (QA/R-5 ELEMENT A9)

The information and documentation methods, data format, and data reporting requirements required by participating laboratories are detailed in *Toxicity Data Submittal Instructions* (**Appendix C**).

All data, reports, and documents submitted to SCCWRP will be stored in electronic format as scans of hard copies and Microsoft Excel files saved in a database that can be made publicly available after completion of the project. Similar to the policy used for the historical data collected, the identity of individual laboratories will be kept anonymous in the database and proprietary information will not be shared with the public.

The SCCWRP Co-Principal Investigator/Project Manager will also ensure that each participating laboratory and SCCWRP personnel participating in this ILS are provided with any updated versions of this QAPP to ensure that ILS protocols and data reporting requirements are consistently followed.

## 2 DATA GENERATION AND ACQUISITION

### 2.1 EXPERIMENTAL DESIGN (QA/R-5 ELEMENT B1)

This section will refer to SCCWRP's *Baseline Testing for Ceriodaphnia dubia Toxicity Testing Laboratory Standardization – Study Plan and Logistics* (**Appendix D**) to describe the approach for ILS testing that was summarized in Section A6 – Project Description.

## 2.2 SAMPLING METHODS (QA/R-5 ELEMENT B2)

This section describes the methods for preparing the samples, dilution waters, and dilution series; subsampling for chemical analysis; splitting and distributing samples.

### 2.2.1 Production of Waters for Testing and Sampling

Two types of water will be produced for this testing. Water made using the EPA moderately hard (EPAMH, sample 1) formula will be used as one of the unspiked dilution water recipes. The EPA MH water will be produced by a contracted laboratory. The water will be made using the formula in the EPA manual without addition of selenium. SCCWRP will pick up the water for each testing round in a 20-gallon plastic carboy. A carboy will be delivered to the laboratory ahead of each round of testing. The total volume needed will be 20 gallons for each round. Once the water is transported to SCCWRP it will be allowed to equilibrate at 4° C in the walk-in fridge for 2 days with aeration using clean air stones, new tubing and portable air pumps before being transferred to cubitainers for the individual laboratories.

A second dilution water recipe will be prepared at SCCWRP using the “Perrier” method. This method will be used to produce sample 2A (unspiked) and to prepare a dilution series using sodium chloride as the spiking agent (samples 2B-F). The Perrier water will be made in a large plastic tank with a 20:80 mix of Perrier to MilliQ water. The water will be stirred on a large vortex mixer and aerated overnight. The water will then be equilibrated in the 4° C walk-in fridge for an additional day before being spiked with sodium chloride. The total needed for all samples using this recipe will be about 100 gallons. This will require approximately 20 gallons of Perrier.

### 2.2.2 Spiking

The most concentrated sample of the dilution series (Sample 2B; 2000 mg/L NaCl) will be created by dissolving NaCl (VWR Life Science, High Purity Grade, CAS 7647-14-5) into the “Perrier” water. The spiking will be done in 20-gallon carboys fitted with a spigot. Four additional concentrations will be produced (samples 2C-F; 1000, 500, 250, 125 mg/L NaCl) using a 1:1 dilution (i.e., 50%). The samples will be mixed for an hour on a large vortex mixer and then allowed to equilibrate in the 4°C ± 2 in the walk-in fridge for ~24 hours before being transferred to cubitainers for the individual laboratories.

For sample 3, 7.00 g of NaCl (same source as above) will be weighed and placed in 100 ml HDPE containers. Each laboratory will receive one container per round. Laboratories will be instructed to dilute the supplied NaCl in 3.5 L with their control dilution water and perform a 50% dilution series to generate a total of 5 dilutions (2000, 1000, 500, 250, and 125 mg/L).

## 2.3 SAMPLE HANDLING AND CUSTODY (QA/R-5 ELEMENT B3)

This section describes sample handling procedures and considerations.

Samples will be shipped to each laboratory (see Appendix A) according to the schedule agreed upon with the participating laboratories. Samples will be shipped on wet ice using priority overnight courier service (i.e., OnTrac or FedEx Priority Overnight). The shipments will also include chain-of-custody (COC) forms completed by SCCWRP, which will reference the study plan and testing instructions. SCCWRP will notify the laboratories via email once the samples are in transit and provide a tracking number.

The cubitainers must be maintained at  $4^{\circ}\text{C} \pm 2$  prior to use and may be held for up to 48 hours after being collected from the parent sample before first use.

If a sample is not delivered to a laboratory by 2 p.m. on the expected arrival date or if the sample has spilled during shipment, the laboratory must contact SCCWRP promptly. SCCWRP will ship new samples that same day. These replacement samples must be used to begin tests within 24 hours of receipt (i.e., within 48 hours of when the initial samples were collected from the original bulk sample).

A laboratory will be given an opportunity to retest a sample if acceptability criteria are not met. However, both the failed test and retest data must be submitted to SCCWRP. A laboratory planning to retest must contact SCCWRP within 24 h of knowing that a test failed acceptability criterion. Laboratories are encouraged to retest with remaining sample. However, arrangements could be made to retest with the archived sample provided by SCCWRP, if possible.

## 2.4 ANALYTICAL METHODS (QA/R-5 ELEMENT B4)

### 2.4.1 Toxicity Tests

*C. dubia* toxicity testing will be performed according to EPA (2002b, 2016) guidance for test method 1002.0. Specific test conditions are described in Table 2. Each laboratory will follow its own specific SOPs and testing protocols for this test with the exceptions described in Section 1.6.2 of Appendix D. Note if any SOPs have been updated since July 2021, laboratories shall provide the updated documents to SCCWRP. Study-specific procedures are intended only to collect additional data (e.g., extending the test through day-8 regardless of when the typical test termination criterion is met) and should not result in any changes to the normal test implementation. Labs must count and report the live young and may report any dead young separately as well. Any dead young that are noted should be reported in the comments field of data submission template. Labs should also document when split broods are detected and report their observations in the comments field of data submission template. An animal that dies before producing young, if not identified as a male, will be included in the analysis with zero (0) to



represent the number of neonates. Thus, the resulting number of live neonates per female combines the survival and reproduction endpoints.

Toxicity test results will be described by SCCWRP as the percent survival and neonates per female. Concentration response relationships may also be described for dilution series by determining the no observed effect concentration (NOEC), the lowest observed effect concentration (LOEC), and a point estimates to determine the IC25 and other relevant percent effect levels as described by EPA (2002b, Appendix M). The type of concentration response, such as those describe by EPA (2000b), will also be identified, and reported.

## 2.4.2 Analyses of Water

Water quality measurements (i.e., conductivity, DO, pH, and temperature) by each toxicity testing laboratory will be conducted according to each laboratory's SOPs and conform with the quality criteria described in Table 3. Water samples submitted to an analytical laboratory will be analyzed according to the methods, target MDLs, and target RLs listed in **Table 4**.

The method detection limit (MDL) is the minimum concentration of a substance that can be measured and reported with 99 percent confidence that the analyte concentration is greater than zero, as determined from analysis of a sample in a given matrix containing the analyte. The analytical laboratory will calculate and report an MDL for each analyte of interest in each matrix (i.e., water) before analyzing samples. The laboratory will calculate the MDLs statistically, based on instrument performance, at least once annually for each analytical method. Ideally, the reporting limits (RLs) will be at least two times the calculated MDL to assure that the quantification of detected compounds is valid. The laboratory will report any results between the MDL and the RL as estimated and qualified by a J-flag. Quantification based on extrapolation will not be accepted. If samples are outside the calibration range, they must be diluted or concentrated, as necessary, and reanalyzed.

## 2.5 QUALITY CONTROL (QA/R-5 ELEMENT B5)

This section describes QC sample requirements and their frequency of collection. Error! Reference source not found. presents laboratory quality control samples for *C. dubia* toxicity testing. **Table 5** presents quality control samples for chemical analyses.

**Table 5. Quality control samples for baseline interlaboratory *C. dubia* toxicity tests.**

Quality Control Sample	Description	Required	Sampling Frequency / Rationale
Negative Control	Laboratory controls will be prepared by each toxicity laboratory to (1) assess acceptable performance of the organism in the absence of a contaminant and (2) compare to organism performance in test samples.	Yes	Concurrent testing with each batch at each toxicity laboratory is required by the EPA (2002b) method.
Positive Control	A reference toxicant test is a toxicity test with a dilution series with a consistent toxicant (e.g., NaCl) in standard dilution water. The test is repeated over time to assess the consistency or organism responses.	No	Participating laboratories are using different chemicals at different concentrations. Instead of running their reference toxicant, all laboratories will provide their control charts for the last 12 months (i.e., Oct 2021 to Sept 2022).

**Table 5. Quality control (QC) samples for chemical analyses.**

QC Sample	Description	Required	Sampling Frequency / Rationale
Duplicates	Samples collected to monitor the precision of sampling.	No	A duplicate for chemical analyses will be collected and analyzed at a frequency of 5% of the total project sample count.
Equipment Blank	Samples of deionized water passed through and/or over the surface of decontaminated sampling equipment to evaluate the potential for sample contamination from equipment.	No	No equipment is needed. Samples for chemical analysis should be decanted directly in the sample containers.
Laboratory Blank	To assess potential ambient sources of contamination associated with decontamination and sampling procedures	Yes	Laboratory blanks will be analyzed by the analytical lab where appropriate.
Matrix Spike/ Matrix Spike Duplicate	Spike recoveries are used to evaluate potential matrix interferences, as well as accuracy. The duplicate spike results (MS and MSD) are compared to evaluate precision.	Yes	Analytical precision and accuracy will be evaluated using laboratory derived matrix-specific MS/MSDs unless there is a potential for matrix interference in chemical analyses.

### 2.5.1 Toxicity Test

A laboratory control for toxicity testing is a negative control that describes *C. dubia* performance absent chemical stressors. This is required with each toxicity test to evaluate test acceptability criteria, as a basis for determining the relative percent effect in a sample, and for determining if there are statistically significant effects in a sample. The control also provides information on

stock organism health and the normal variability in survival and reproduction of those stock organisms (SWRCB 2020b).

Each test will be run for eight days with all neonate production reported by the laboratories. After QA evaluation of the data submitted, SCCWRP will determine the production of three broods for each test in a consistent manner. Total neonate production for each female will then be calculated from these three broods and used in subsequent analysis.

Variability of toxicity test results will be evaluated using repeated analyses of the same samples within a single laboratory to assess intra-laboratory precision (i.e., repeatability). This will be accomplished by repeated testing among multiple batches. Likewise, repeated analyses of the same samples among laboratories will be compared to assess inter-laboratory precision (i.e., reproducibility). The mean and CV among replicate analyses will be used to describe variability.

### **2.5.2 Chemical Analysis**

Duplicate samples, blanks, and matrix spike/matrix spike duplicates (MS/MSDs) will be submitted for chemical analyses as indicated in **Table 6**. SCCWRP will randomly identify a participating toxicity laboratory or laboratories to collect these samples. Equipment blanks are not required because no equipment will be needed to collect samples.

## **2.6 INSTRUMENT/EQUIPMENT TESTING, INSPECTION, AND MAINTENANCE (QA/R-5 ELEMENT B6)**

Analytical instrument testing, inspection, maintenance, setup, and calibration will be conducted by the laboratories in accordance with the requirements identified in the laboratory's SOPs and manufacturer instructions. Each analytical method provides protocols for proper instrument setup, tuning, and critical operating parameters that will be followed. Instrument maintenance and repair will be documented by each laboratory in maintenance logs or record books to be available upon request.

## **2.7 INSTRUMENT/EQUIPMENT CALIBRATION AND FREQUENCY (QA/R-5 ELEMENT B7)**

Laboratory instruments will be properly calibrated, and the calibration will be verified with appropriate check standards and calibration blanks for each parameter before beginning each analysis. Instrument calibration procedures and schedules will conform to analytical protocol requirements and descriptions in the SOP of each participating laboratory.

Calibration standards will be obtained from a commercial vendor traceable back to the National Institute of Standards and Technology (NIST) where available and appropriate. Stock standards will be used to establish intermediate standards and calibration standards. Special attention will

be given to expiration dates, proper labeling, proper refrigeration and storage, and prevention of contamination. Documentation relating to the receipt, mixing, and use of standards will be recorded in a laboratory logbook. All analytical laboratory calibration and spiking standards will be confirmed against standards from another source, or as specified in the testing method and laboratory SOP. This data may be requested from the laboratories at a later date.

## 2.8 INSPECTION/ACCEPTANCE OF SUPPLIES AND CONSUMABLES (QA/R-5 ELEMENT B8)

The quality of supplies and consumables used during the baseline ILS can affect data quality. All equipment that comes into contact with samples must be sufficiently clean and absent residues to prevent detectable contamination. Analyte concentrations must also be accurate in all standards used for calibration and quality control purposes.

Supplies and consumables needed to complete this baseline ILS include the following:

- Sufficiently large carboys are needed to hold the entire volume of each sample for all testing and analyses so that each bulk sample for each batch can be prepared in a single container.
- Deionized water or Milli-Q water® and reagent grade chemicals will be used to prepare synthetic dilution waters as described by EPA (2002b, 2016). These chemicals include NaHCO<sub>3</sub>, CaSO<sub>4</sub>\*2H<sub>2</sub>O, MgSO<sub>4</sub>, KCl, NaCl, and Perrier® water (if demineralized water is required).
- New one (1) gallon cubitainers will be used to transport split samples to participating laboratories for toxicity testing.
- Coolers and double-bagged ice will be used to transport samples to each participating laboratory.
- *C. dubia* for toxicity testing will be provided from in-house cultures by each participating laboratory. Cultures producing fewer than 20 young/female or with adult mortality greater than 20% mortality will not be used (EPA 2002b, 2016).
- Control water and food recipes prepared according to each participating laboratory's SOPs, *C. dubia* food prepared according to each participating laboratory's SOPs will be provided by each participating laboratory.
- New sample containers for chemical analysis will be provided by SCCWRP to each participating laboratory.
- Laboratory water and cleaning solutions used for decontamination will adhere to standard SOPs for each laboratory and will be documented.

## **2.9 NON-DIRECT MEASUREMENTS (QA/R-5 ELEMENT B9)**

Existing chemical and toxicity data from prior toxicity tests will be used for this study. Reference toxicity testing to determine the running mean and quality limits (i.e., 2 standard deviations from the running mean) from the last 12 months will be used to describe the normal range of *C. dubia* responses to a consistent stressor at each laboratory. Any existing data will be reviewed for quality assurance and acceptability prior to use in the SWB *C. dubia* study.

## **2.10 DATA MANAGEMENT (QA/R-5 ELEMENT B10)**

Data for this baseline ILS will be generated at analytical laboratories. Data developed by or provided to SCCWRP by participating laboratories will be maintained by SCCWRP electronically in a server location dedicated to the SWB *C. dubia* study. Files will include complete data reports provided by analytical and toxicity testing laboratories, copies of any sample preparation logs, COCs, toxicity testing bench sheets, photos, and communications among laboratory participants. Data files will be routinely backed-up and a copy of current data will be stored securely at a separate location maintained by SCCWRP.

# **3 ASSESSMENT AND OVERSIGHT**

## **3.1 ASSESSMENT AND RESPONSE ACTIONS (QA/R-5 ELEMENT C1)**

Laboratories will promptly report any deviations from this QAPP to SCCWRP. The need for corrective actions will then be determined by the SCCWRP Technical Coordinator, SCCWRP Project Manager, SCCWRP QA Officer, Laboratory Coordinator, and Laboratory QA Manager by evaluating the root cause of any deviations from this QAPP and determining the potential effects of the deviation on data quality. Corrective actions will provide specific guidance and/or requirements to remedy the cause of such deviations. The SCCWRP QA Officer will provide a summary of the deviation, root cause, and any needed corrective action to the SCCWRP project team before additional sampling is conducted. Corrective actions could include resampling and/or re-analyzing samples if, for example, hold times have been exceeded, if data quality indicators are not met, or if toxicity tests do not meet test acceptability criteria.

## **3.2 REPORTING (QA/R-5 ELEMENT C2)**

Data developed from this baseline ILS will be summarized and reported to the SAC and ESP and made available to the public. Initial toxicity and chemical analysis results from interlaboratory

testing will be summarized by SCCWRP and reported to the ESP and SAC within 60 business days following the complete submission of the baseline testing data by the laboratories.

SCCWRP will provide a description of the tested samples, the range of toxicity test results (e.g., live neonates per female, percent difference from controls, IC25, etc.), and the range of test variables (e.g., water quality characteristics, dilution water type, and specific 8-hour age range of neonates within the <24-hour old test requirement) will be summarized and compared within and among laboratories. SCCWRP will describe factors associated with inter- and intra-laboratory test variability that could be the subject of targeted experiments to standardize select test parameters to minimize inter- and intra-laboratory variability, if any are identified. Any laboratory tests that do not meet test acceptability criteria or results that do not meet MQOs and are qualified or rejected will also be presented. Deviations from the QAPP by any laboratory, lessons learned, and any corrective actions will also be discussed in the report.

## **4 DATA REVIEW AND USABILITY**

This section describes the data verification and validation procedures that occur after the data collection or generation phase of the project to ensure that data are sufficient to meet the project goals and objectives by complying with specified criteria.

### **4.1 DATA REVIEW, VERIFICATION, AND VALIDATION (QA/R-5 ELEMENT D1)**

A review to verify and validate data quality will be performed after each sampling event. These evaluations will identify any deviations from the QAPP and confirm that MQOs were met (e.g., COC forms were complete, sampling hold times were met, RLs were met). Corrective actions may be implemented if data requirements are not met.

### **4.2 VERIFICATION AND VALIDATION METHODS (QA/R-5 ELEMENT D3)**

Data verification, confirming that the correct values are recorded and are transcribed accurately into EDFs, will be performed initially by each laboratory prior to transmitting any data packages to SCCWRP. Secondary data verification will be conducted by SCCWRP by spot checking data submittals at a rate of 5% of the toxicity results (i.e., 1 of every 20 electronic data records will be confirmed by checking the reported value with the toxicity bench sheet) to confirm that bench sheet data are correctly reflected in the electronic files or summary tables provided by testing laboratories. These checks will also confirm the following:

- Required analyses were conducted

- Expected methods were used
- Target RLs and MDLs were met
- Batch and sample identifications were reported correctly
- Units were reported correctly
- Quality control evaluation calculations were complete
- COCs were complete
- Data were reported according to the required format

Any data inconsistencies or questions identified by SCCWRP will be discussed with the laboratory and either a corrected data file or documentation of the explanation will be recorded in the project file.

Data validation will be performed by SCCWRP following the receipt of each batch of toxicity test data to ensure that the data developed meet data quality requirements and are adequate to meet the study goals. Validation will evaluate the following:

- Hold times were met
- Detections in lab control samples
- Exceedances of MQOs
- Toxicity test acceptability criteria were met
- Toxicity test requirements specific to this ILS (e.g., test set-up, duration, and data requirements) were met

Any toxicity test results where the controls do not meet test acceptability criteria defined in Table 2 will be rejected for use as valid tests. Although, these test results may still be used to identify test parameters that contribute to inter- and intra-laboratory variability and may be included in statistical analyses to meet the study goals. Likewise, toxicity data where all test requirements were not met shall be qualified but may still be included in analyses to identify test parameters that contribute to inter- and intra-laboratory variability. Causes of any failures to meet test acceptability criteria or test requirements will be reported to SCCWRP by the laboratory. The SCCWRP Technical Coordinator, SCCWRP Project Manager, SCCWRP QA Officer, Laboratory Coordinator, and Laboratory QA Manager will determine if and how qualified results may be used and SCCWRP will document any corrective actions that were implemented to prevent failures from recurring.

### 4.3 RECONCILIATION WITH USER REQUIREMENTS (QA/R-5 ELEMENT D3)

Data developed from this ILS will be evaluated to identify factors contributing to inter- and intra-laboratory test variability. Such factors, if identified, may inform targeted experiments to standardize test parameters with the intent of minimizing inter- and intra-laboratory variability as part of Task 3 of the SWB *C. dubia* study.



## 5 REFERENCES

- American Society for Testing and Materials (ASTM). 1999. Standard practice for conducting an interlaboratory study to determine precision of a test method, E691-20. In: Annual Book of ASTM Standards, Vol. 14.05. Philadelphia, PA.
- Central Valley Clean Water Agency (CVCWA). 2018. Toxicity Special Study: Phase I Study Report. Prepared by Larry Walker Associates, Robertson-Bryan, Inc., and Pacific EcoRisk. December 2018.  
[https://www.waterboards.ca.gov/water\\_issues/programs/state\\_implementation\\_policy/docs/comments20181221/debbie\\_webster.pdf](https://www.waterboards.ca.gov/water_issues/programs/state_implementation_policy/docs/comments20181221/debbie_webster.pdf).
- Schiff, K.C. and D. Greenstein. 2016. Stormwater Monitoring Coalition Toxicity Testing Laboratory Guidance Document. Southern California Coastal Water Research Project (SCCWRP), Technical Report 956. Available:  
[https://ftp.sccwrp.org/pub/download/DOCUMENTS/TechnicalReports/956\\_StrmWtrMonitCoalitToxTestingLabGuid.pdf](https://ftp.sccwrp.org/pub/download/DOCUMENTS/TechnicalReports/956_StrmWtrMonitCoalitToxTestingLabGuid.pdf)
- Southern California Coastal Water Research Project (SCCWRP). 2022. Stakeholder Meeting Presentation. February 23, 2022. Available: <https://www.sccwrp.org/about/research-areas/additional-research-areas/ceriodaphnia-toxicity-testing-quality-assurance/>
- Southern California Coastal Water Research Project (SCCWRP). 2021a. Stakeholder Meeting Presentation. January 20, 2021 Available: <https://www.sccwrp.org/about/research-areas/additional-research-areas/ceriodaphnia-toxicity-testing-quality-assurance/>
- Southern California Coastal Water Research Project (SCCWRP). 2021b. Conceptual Workplan: Quality Assurance Evaluation of the *Ceriodaphnia dubia* Reproduction Test. Draft. 5 m/Prepared for the State Water Resources Control-087bnm/Board. Agreement #19-278-0780. May 26, 2021. Available:  
<https://ftp.sccwrp.org/pub/download/DOCUMENTS/CeriodaphniaQA/May2021DraftWorkplan.pdf>
- Southern California Coastal Water Research Project (SCCWRP). 2020. October 30, 2020 – Development of Quality Assurance Recommendations for the *Ceriodaphnia dubia* Toxicity Test – Excerpt from contract between the State Water Resources Control Board and the Southern California Coastal Water Research Project Authority (Agreement Number 19-078-270). Available:  
[https://www.waterboards.ca.gov/water\\_issues/programs/state\\_implementation\\_policy/txass\\_cntrl.html](https://www.waterboards.ca.gov/water_issues/programs/state_implementation_policy/txass_cntrl.html)
- State Water Resources Control Board (SWB). 2020a. Development of Quality Assurance Recommendations for the *Ceriodaphnia dubia* Toxicity Test – Excerpt from contract between the State Water Resources Control Board and the Southern California Coastal Water Research Project Authority (Agreement Number 19-078-270). October 30, 2020.

Available:

[https://www.waterboards.ca.gov/water\\_issues/programs/state\\_implementation\\_policy/tx\\_ass\\_cntrl.html](https://www.waterboards.ca.gov/water_issues/programs/state_implementation_policy/tx_ass_cntrl.html)

State Water Resources Control Board (SWB). 2020b. Staff Report, Including Substitutes Environmental Documentation for the Establishment of the Water Quality Control Plan for Inland Surface Waters, Enclosed Bays, and Estuaries of California; and Toxicity Provisions. Available:

[https://www.waterboards.ca.gov/water\\_issues/programs/state\\_implementation\\_policy/tx\\_ass\\_cntrl.html](https://www.waterboards.ca.gov/water_issues/programs/state_implementation_policy/tx_ass_cntrl.html)

State Water Resources Control Board (SWB). 2020c. Surface Water Ambient Monitoring Program Quality Assurance Program Plan. State Water Resources Control Board. Version 1.3. May 2020.

[https://www.waterboards.ca.gov/water\\_issues/programs/swamp/quality\\_assurance.html](https://www.waterboards.ca.gov/water_issues/programs/swamp/quality_assurance.html)

United States Environmental Protection Agency (EPA). 2001a. Final report: interlaboratory variability study of EPA short-term chronic and acute whole effluent toxicity test methods, Vol. 1. Office of Water, U.S. Environmental Protection Agency, Washington, D.C. 20460. EPA/821/B-01/004.

EPA. 2001b. Final report: interlaboratory variability study of EPA short-term chronic and acute whole effluent toxicity test methods, Vol. 2: Appendix. Office of Water, U.S. Environmental Protection Agency, Washington, D.C. 20460.

EPA. 2002a. Guidance for Quality Assurance Project Plans (EPA QA/G-5), EPA/240/R-02/009. Washington, D.C. Available at: <https://www.epa.gov/sites/production/files/2015-06/documents/g5-final.pdf>

EPA. 2002b. Short-term Methods for Estimating the Chronic Toxicity of Effluents and Receiving Waters to Freshwater Organisms. Fourth Edition. Office of Water. Washington, DC. EPA-821-R-02-013. October. <https://www.epa.gov/cwa-methods/chronic-toxicity-freshwater-wet-methods>

EPA. 2006. Guidance on Systematic Planning using the Data Quality Objectives Process (EPA QA/G-4). Washington, D.C. <https://www.epa.gov/sites/production/files/2015-06/documents/g4-final.pdf>

EPA. 2016. Whole Effluent Toxicity Methods Errata Sheet. 28 p. Office of Water. Washington, DC. EPA 821-R-02012-ES. December 2016 <https://www.epa.gov/cwa-methods>.

**Appendix A. Contact information for laboratories participating in inter-laboratory *C. dubia* chronic toxicity testing as part of the California State Water Resources Control Board *C. dubia* study.**

<b>ID#</b>	<b>Laboratory</b>	<b>Contact Name</b>	<b>Contact information</b>	<b>Shipping Address</b>
1.	Aquatic Bioassay & Consulting Labs, Inc	Joe Freas	805-643-5621 x18 joe@aquaticbioassay.com	29 N Olive St., Ventura, CA 93001
2.	Aquatic Testing Laboratory	Joe LeMay	805 650-0546 jlemay12@pacbell.net	4350 Transport Street, Unit 107 Ventura, CA 93003
3.	AquaScience	Kimberly Miller	530-753-5456 Kimberley@aquascience.com	630 Cantrill Dr., Davis, CA 95618
4.	Enthalpy	Peter Arth	858-587-7333 ext. 214 Peter.arth@enthalpy.com	4340 Vandever Avenue, San Diego, CA 92120
5*.	GEI	Natalie Love	303-264-1070 Nlove@geiconsultants.com	4601 DTC Boulevard, Suite 900, Denver, CO, 80237
6*.	Inland Empire	Sushmitha Reddy	909-993-1813 Sreddy@ieua.org	IEUA Water Quality Laboratory, Building C, 6075 Kimbal Ave., Chino, CA 91708
7.	McC Campbell	Drew Gantner	925-252-9262 Drew.gantner@mcccampbell.com	1534 Willow Pass Road Pittsburg, CA 94565-1701
8.	MBC Applied Environmental Sciences	Sonja Beck	714-850-4830 x225 Sbeck@mbcaquatic.com	3000 Redhill Ave., Costa Mesa CA, 92626
9.	Pacific Ecorisk	Stephen Clark	707-207-7760 Slclark@pacificecorisk.com	2250 Cordelia Road. Fairfield, CA 94534
10.	Sanitation Districts of Los Angeles County	Josh Westfall	562-908-4288 x2815 Jwestfall@lacsdc.org	San Jose Creek Biology Lab. 1965 Workman Mill Rd. Whittier, CA 90601
11*.	TetraTech	Marcus Bowersox	410-902-3142 Marcus.Bowersox@tetrattech.com	10711 Red Run Blvd., Suite 105, Owings Mills, MD 21117
12.	Wood	Steve Carlson	858-299-5368 Steve.carlson@woodplc.com	4905 Morena Blvd. Ste. 1304, San Diego, CA 92117
Note: * Laboratory is not currently ELAP accredited				

**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

Step	DQO Guidance of Purpose and Outputs of Step	<i>C. dubia</i> Interlaboratory Study
1. State the Problem	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>Define the problem that necessitates the study; identify the planning team, examine budget, and schedule.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>A concise description of the problem.</li> <li>A list of the planning team members and identification of the decision makers.</li> <li>A description of overall approach for assessment.</li> <li>A summary of available resources, constraints, and relevant deadlines for the study.</li> </ul>	<p><b>Problem</b></p> <ul style="list-style-type: none"> <li>Variability in the test outcome can occur when split samples are submitted to two or more laboratories for chronic toxicity testing with <i>C. dubia</i> so that effluent or blank spit samples have been reported to cause a significant effect at one laboratory and no significant effect at another laboratory.</li> </ul> <p><b>Decision-makers and team members</b></p> <ul style="list-style-type: none"> <li>The SWB is funding and managing a study titled Development of Quality Assurance Recommendations for the Ceriodaphnia Toxicity Test (SWB <i>C. dubia</i> Study) as part of implementing Resolution No. 2020-0044 to establish the Water Quality Control Plan for Inland Surface Waters, Enclosed Bays, and Estuaries of California and Adopting Toxicity Provisions (Toxicity Provisions).</li> <li>SCCWRP is implementing this work on behalf of the SWB in coordination with commercial and municipal laboratories.</li> <li>A Stakeholder Advisory Committee (SAC) and Expert Science Panel (ESP) provide review and input on the study design and reporting.</li> </ul> <p><b>Overall approach</b></p> <ul style="list-style-type: none"> <li>Baseline interlaboratory testing with split samples (the interlaboratory study or ILS) is being conducted as part of the overall SWB <i>C. dubia</i> Study to assess factors contributing inter- and intra-laboratory variability.</li> </ul> <p><b>Resources</b></p> <ul style="list-style-type: none"> <li>The SWB is funding and managing the overall SWB <i>C. dubia</i> study.</li> <li>ILS testing performed by participating laboratories is being funded by CASA.</li> </ul> <p><b>Constraints</b></p> <ul style="list-style-type: none"> <li>The capacity of participating laboratories to perform ILS testing may be limited.</li> <li>Scheduling tests among participating laboratories will need to be coordinated.</li> <li>Available funding may limit the number of tests that can be performed.</li> </ul> <p><b>Deadlines</b></p>

**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

Step	DQO Guidance of Purpose and Outputs of Step	<i>C. dubia</i> Interlaboratory Study
		<ul style="list-style-type: none"> <li>The baseline ILS results were to be reported in a technical memorandum by January 1, 2023.</li> <li>A Final Recommendations Report for the SWB <i>C. dubia</i> study is to be completed by March 30, 2023.</li> </ul>
<p>2. Identify the Goals of the Study</p>	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>State how environmental data will be used in meeting objectives and solving the problem, identify study questions, define alternative outcomes.</li> </ul> <p><b>Approach</b></p> <ul style="list-style-type: none"> <li>Identify the key question that the study attempts to address and alternative outcomes or actions that may be taken, depending on the answer to the key study question; and develop decision statements.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>A statement of the decision that must be resolved using data in order to address or solve the problem.</li> <li>A list of possible actions or outcomes that would result from each resolution of the decision statement.</li> </ul>	<p><b>Overall goal</b></p> <ul style="list-style-type: none"> <li>Identify potential sources of variability within and among laboratories and make recommendations to the SWB regarding how the chronic <i>C. dubia</i> test method could be conducted to reduce sources of variability in the results (SWB 2020a, SCCWRP 2021b).</li> </ul> <p><b>Key questions</b></p> <ul style="list-style-type: none"> <li>Which lab practice(s) should be standardized to reduce intra- and inter-lab variability?</li> <li>Does standardizing lab practices improve consistency and comparability in <i>C. dubia</i> test results within and among laboratories?</li> </ul> <p><b>Possible outcomes</b></p> <ul style="list-style-type: none"> <li>One or more controllable variables are identified that are associated with inter- and/or intra-laboratory variability and recommendations can be made to the SWB to reduce this source(s) of variability.</li> <li>One or more variables are identified that are associated with inter- and/or intra-laboratory variability but they are not controllable and recommendations cannot be made to the SWB to reduce this source(s) of variability.</li> <li>No consistent source(s) of variability associated with the method can be identified and recommendations cannot be made to the SWB to reduce source(s) of variability in test results.</li> </ul>
<p>3. Identify Information Inputs</p>	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>Identify the data and information needed to answer study questions.</li> </ul> <p><b>Activities</b></p>	<p><b>Informational inputs</b></p> <ul style="list-style-type: none"> <li><i>C. dubia</i> chronic toxicity testing with splits of positive and negative control samples among multiple laboratories in California.</li> </ul>

**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

Step	DQO Guidance of Purpose and Outputs of Step	<i>C. dubia</i> Interlaboratory Study
	<ul style="list-style-type: none"> <li>• Identify the types and sources of information needed to resolve decisions or produce estimates.</li> <li>• Identify the sources for each item of information identified.</li> <li>• Select appropriate sampling and analysis methods to provide the necessary data.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>• A list of informational inputs (including sources and potential action levels) needed to resolve the decision.</li> <li>• The list of environmental variables or characteristics that will be measured.</li> </ul>	<ul style="list-style-type: none"> <li>• Repeated analyses within participating laboratories (i.e., 2 or 3 batches of baseline testing).</li> <li>• Water quality data and other testing or culture data collected during split-sample testing.</li> </ul> <p><b>Variables/characteristics to be measured</b></p> <ul style="list-style-type: none"> <li>• Chemical analysis of split samples to document consistency or test materials (alkalinity, conductivity, hardness, major ion concentrations, pH).</li> <li>• Water quality analysis of tested samples to document conditions during testing (alkalinity, conductivity, dissolved oxygen, hardness, pH, and temperature).</li> <li>• Laboratories will report: daily survival and neonates per female; number of males, unhealthy and dead adults, and any dead neonates in the brood board; specific beginning and end time window for age of neonates at test initiation; dilution water recipe and food recipe; and, light intensity and air temperature within the testing area at the time of the testing.</li> </ul>
<p>4. Define the Boundaries of the Study</p>	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>• Specify the target population and characteristics of interest; define the spatial and temporal limits, and scale of inference.</li> </ul> <p><b>Activities</b></p> <ul style="list-style-type: none"> <li>• Define the target populations/parameters (physical, chemical, and biological), as applicable.</li> <li>• Specify the spatial and temporal boundaries.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>• Characteristics that define the domain of the assessment.</li> <li>• A detailed description of the spatial and temporal boundaries of the characterization.</li> </ul>	<p><b>Study domain</b></p> <ul style="list-style-type: none"> <li>• Participating laboratories (n= 4) are among those currently or recently accredited by ELAP to perform the <i>C. dubia</i> chronic toxicity test in California.</li> <li>• Baseline toxicity testing will be conducted consistent with EPA (2002) guidance and according to each laboratory's SOP.</li> </ul> <p><b>Spatial domain</b></p> <ul style="list-style-type: none"> <li>• Participating laboratories (10 out of 12) are located throughout California and two laboratories are out of state in Colorado and Maryland.</li> </ul> <p><b>Temporal domain</b></p> <ul style="list-style-type: none"> <li>• Testing is to be performed in 2022 over a 8-week period.</li> </ul>

**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

Step	DQO Guidance of Purpose and Outputs of Step	<i>C. dubia</i> Interlaboratory Study
5. Develop the Analytical Approach	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>Define the parameter of interest, specify the type of inference, and develop the logic for drawing conclusions from findings.</li> </ul> <p><b>Activities</b></p> <ul style="list-style-type: none"> <li>Specify the parameters (e.g., mean, median, percentile) that are considered important to characterize the population of interest.</li> <li>Develop decision rules (“if ... then ...”) for the parameters of interest.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>Identify the most relevant population parameters for estimation.</li> <li>Decision rules for interpretation of results.</li> </ul>	<p><b>Parameters of interest</b></p> <ul style="list-style-type: none"> <li>Mean neonates per female in tested samples and controls as well as calculated IC25 and IC50 values will be compared to determine if they are significantly different within and among laboratories.</li> <li>Culturing information (e.g., culture water type and water quality parameters, food quality, culture age) and ancillary testing information (e.g., specific age of test organisms within an 8-hour window, specific time at test termination, etc.).</li> <li>Test data/factors will be evaluated to determine if it is associated with variability in <i>C. dubia</i> survival or reproduction within or among laboratories.</li> </ul> <p><b>Decision rules</b></p> <ul style="list-style-type: none"> <li>Standard test method statistics (EPA 2002b, 2016) may be used to evaluate statistical differences between samples and control endpoints.</li> <li>Percent effect may be calculated for the sample test treatment relative to the laboratory’s standard control water.</li> <li>Inter- and intra-laboratory variability will be evaluated by direct comparisons of the IC values and CVs among repeated analyses.</li> <li>ASTM (1999) h and k statistics will be used to evaluate data consistency and identify potential outliers among laboratories and within a laboratory, respectively.</li> <li>Multi-variate and other statistical analyses will be used to assess test factors contributing to inter- and intra-laboratory variability.</li> </ul>
6. Specify Performance or Acceptance Criteria	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>Develop the performance criteria for new data being collected or acceptable criteria for existing data being considered for use.</li> </ul> <p><b>Activities</b></p>	<p><b>Performance criteria</b></p> <ul style="list-style-type: none"> <li>Toxicity test results must meet test acceptability criteria described by the test method (EPA 2002b, 2016) to be considered valid for comparisons among valid test results.</li> <li>Tests that do not meet test acceptability criteria and all test requirements will also be assessed to determine factors contributing to the failure, which may inform the analysis of variables causing or contributing to inter- and intra-laboratory variability.</li> </ul>

**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

Step	DQO Guidance of Purpose and Outputs of Step	<i>C. dubia</i> Interlaboratory Study
	<ul style="list-style-type: none"> <li>Specify the performance metrics and acceptable levels of uncertainty.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>The primary output is a set of acceptance criteria that collected data should achieve in order to minimize the possibility of failing to keep uncertainty within acceptable limits.</li> </ul>	<ul style="list-style-type: none"> <li>Water quality (e.g., dissolved oxygen, pH, etc.) and sample chemistry (e.g., major ion concentration) data must meet measurement quality objectives defined in Section A7 of this QAPP.</li> </ul>
<p>7. Develop the Plan for Obtaining Data</p>	<p><b>Purpose</b></p> <ul style="list-style-type: none"> <li>Select the resource-effective sampling and analysis plan that meets the performance criteria.</li> </ul> <p><b>Activities</b></p> <ul style="list-style-type: none"> <li>Compile information that we will need (generated in Steps 1 through 6) in an acceptable and efficient sampling and analysis design.</li> <li>Identify constraints that will affect the sampling and analysis design.</li> <li>Provide details on the sampling and analysis methods that we will use to generate the data.</li> <li>Prepare a resource-effective information collection plan that will meet our needs and requirements.</li> </ul> <p><b>Outputs from this step</b></p> <ul style="list-style-type: none"> <li>The most resource-effective design for the study that is expected to achieve the DQOs, selected from a group of alternative designs.</li> </ul>	<p><b>Approach</b></p> <ul style="list-style-type: none"> <li>The sampling and analysis approach is described in this QAPP and relies on standardized methods for data collection and analysis.</li> <li>SCCWRP and the ESP are currently discussing the study design and this section will be updated with study design details once determined.</li> </ul>

Notes:

CASA – California Association of Sanitary Agencies  
 ELAP – Environmental Laboratory Accreditation Program  
 ESP – Expert Science Panel  
 QAPP – Quality Assurance Project Plan



**Appendix B. Data quality objectives for the SWB *C. dubia* study interlaboratory testing.**

---

<b>Step</b>	<b>DQO Guidance of Purpose and Outputs of Step</b>	<b><i>C. dubia</i> Interlaboratory Study</b>
-------------	--	--

---

SAC – Stakeholder Advisory Committee  
SCCWRP – Southern California Coastal Watershed Research Project  
SOP – Standard Operating Procedure  
SWB – California State Water Resources Control Board

## Appendix C – Toxicity Data Submittal Instructions

This section provides guidance for the use of the data submission template provided by SCCWRP to submit the results of this ILS. The template is an Excel file with seven tabs corresponding to different data types. Each tab and the fields within are described below. No fields can be left blank. In cases where no data is available, a default value will be provided.

Once filled out, the template must be submitted through the SCCWRP data portal. The portal will perform automatic checks for business rules and verify that the data is complete and meet study specific requirements. Submissions not meeting the rules will be returned automatically with descriptions of the problems with the file. The file can then be corrected and resubmitted. Once the data passes the checking process a receipt will be issued confirming success.

For any questions or issues with data entry and submission, please contact Darrin Greenstein at [darring@sccwrp.org](mailto:darring@sccwrp.org).

### LabTestInfo

Per the laboratories request, this tab will be filled as a survey form. A separate link will be sent to the laboratories. The instructions below pertain to the survey form. LabTestInfo provides information on a wide range of test conditions that are expected to be the same for all tests conducted within a given laboratory. If any of the parameters in this tab/survey change during the course of the study, the laboratory must contact SCCWRP for instructions on how to proceed.

*LabCode*-Enter the code from the list below. During the data submission process the system will anonymize each lab to the code they were assigned for the historic data assessment.

Lab Name	LabCode
Aquatic Bioassay	ABC
AquaScience	AQSC
Aquatic Testing Laboratories	ATL
Enthalpy	ENPY
GEI	GEI
Inland Empire	IEUA

Lab Name	LabCode
LA County Sanitation	LASD
Marine Biological Consultants	MBC
McCambell	MCBL
Pacific EcoRisk	PCR
TetraTech	TRTH
Wood	WOOD

*Contact*-First name and last initial of person to contact with regards to data questions.

*DilutionWaterRecipe*-Code for the method used to make dilution water used for laboratory controls and dilution of sample 3. Choose from the list below.

<b>Description</b>	<b>Code</b>
EPA Moderately Hard	EPAMH
EPA Moderately Hard + Se	EPAMH+Se
EPA Moderately Hard + Vitamins	EPAMH+V
EPA Moderately Hard + Se + Vitamins	EPAMH+Se+V
Mineral Water + DIW. Code is mineral water type followed by the ratio of mineral water to DIW (e.g. Perrier 20:80)	XXXX YY:ZZ
Other water type not mentioned above	Contact SCCWRP

*SourceWater*-Brief description of type of water used as the base for dilution water (e.g., Millipore, RO, Resin).

*YCTSource*-Is the YCT made in-house or purchased? Field can be “In-house”, “ARO”, “ABS”, or other. If other, please contact SCCWRP for a code.

*YCTRecipe*-Brief description of each component used to make the YCT (e.g., Fleschmans Yeast+Trader Joes Wheatgrass+Purina Trout Chow).

*AlgaeSource*-Is the YCT made in-house or purchased? Field can be “In-house”, “ARO”, “ABS”, or other. If other, please contact SCCWRP for a code.

*TestChamberMaterial*-Brief description of material that the test chambers are made from (e.g., Polyethylene, Glass)

*TestChamberVolume*-Volume of test chambers expressed in milliliters.

*VolumeTestSolution*-Volume of sample in the test chambers expressed in milliliters.

*FeedingFrequency*-Number times per day that the test chambers are fed.

*FeedingMethod*-Code for feeding procedure; either directly into the chambers or into solutions before addition to the chambers. “Direct” for the former and “Solution” for the latter.

*YCTConcentration*-Concentration of YCT expressed in milliliters of YCT per milliliter of test solution.

*AlgaeConcentration*-Concentration of algae expressed as number cells per milliliter of test solution.

*Comments*-A place for any comments or clarifications to any of the fields in this tabs.

## **ToxBatchInfo**

This tab contains information for each individual test. The definition of a “batch” is a sample or samples (dilution series) and its associated dilution water control. The expected number of lines in this tab should be equal to the number of samples tested in the study for each laboratory. Some fields are repeated in subsequent tabs to facilitate downstream analyses and merging.

*ToxBatch*-Unique identifier assigned by the laboratory. Each grouping of samples and controls must be assigned its own identifier. The identifier can be alphanumeric.

*LabCode*-Same as for the LabTestInfo.

*SampleHoldTime*-Amount of time between the receipt of the sample and the start of the test in hours or decimal days.

*SampleHoldTimeUnits*-Units for sample hold time. Either “hours” or “days”.

*TestStartDateTime*-Date and time of the start of the test.

*TestEndDateTime*-Date and time of the end of the test.

*Matrix*-Designation for a test as either being for the water samples (W) or (S) for the solid NaCl sample.

*ActualTestDuration*-Duration of the test expressed in hours or decimal days.

*ActualTestDurationUnits*-Units for the test duration. Either “hours” or “days”.

*TestAcceptability*-Code for the overall acceptability of the test. See table below. Contact SCCWRP if it appears that a condition occurred that is not covered by a code in the table.

Acceptability Code	Test Condition Description
A	Test meets all acceptability criteria
C	Reduced number of replicates
D	Control performance criteria not met
E	Sample tested outside of specified holding time
H	Water quality data incomplete
J	Minor deviation in test conditions
K	Incoming sample temperature exceeded limits

*ReferenceBatch*-ToxBatch-Not used for baseline study. Enter NR for all rows.

*LightIntensity*-Measured light intensity at the location of the experiment.

*LightIntensityUnits*-Units for the light intensity measurement. Footcandles or microEinsteins/square meter/second (uE/m<sup>2</sup>/s).

*NeonateMinAge*-Minimum age of neonates used to start the test expressed in hours.

*NeonateMaxAge*-Maximum age of neonates used to start the test expressed in hours.

*SampleStoreCond*-Brief description of how sample was stored. Include temperature and light conditions (e.g., 4°C, dark).

*DilWaterPrepDate*-Date that dilution water was prepared.

*YCTPrepDate*-Date that the YCT used in the test was prepared. For purchased YCT, use the date of arrival for non-frozen or the date thawed for frozen. Note in comments if arrival frozen or not.

*AlgaePrepDate*-Date that the algae was prepared. For purchased algae, use the date of arrival.

*Comments*- A place for any comments or clarifications to any of the fields in this tab. Comments should pertain to the specific ToxBatch on each line.

## **ToxResults**

This is the main data tab where survival and reproduction data are entered. Each line is a combination of replicate, day, and time of observation.

*SampleID*-The ID provided by SCCWRP for each sample. For all laboratory control samples, enter 0000 as the sample ID.

*SampleArrivalDateTime*-The date and time that the sample arrived at the laboratory.

*ToxBatch*-Same value as from the ToxBatchInfo tab.

*LabCode*-Same value as from the LabTestInfo tab.

*Dilution*-Strength of the sample. Full strength sample will be 100. For laboratory control, enter -88.

*Concentration*-Not used for the baseline study. Enter -88 for all sample types.

*ConcentrationUnits*- Not used for the baseline study. Enter NR for all sample types.

*EndPoint*-Test endpoints. Either Repro for neonate production or Survival.

*LabRep*-Replicate number for each sample.

*Day*-Test day.

*Time*-The time of the observation.

*Result*-Data for reproduction or survival. Reproduction data must be entered as counts of neonates at each observation point. For survival, enter 1 for live and 0 for dead adult. After an adult dies, enter 0 for all subsequent days and -88 for reproduction.

*ResultUnits*-Unit of measurement for the result. For both reproduction and survival, the value must be Count.

*QACode*-Similar to AcceptabilityCode in ToxBatchInfo, but at the individual observation level. Choose from the table below. Field cannot be blank. Use comments for any explanation of a code that is not A. Contact SCCWRP for any situation that does not fit one of the descriptions.

Acceptability Code	Test Condition Description
A	Test meets all acceptability criteria
C	Replicate lost or observation not made.
E	Sample tested outside of specified holding time
J	Minor deviation in test conditions

*SampleTypeCode*-Code for one of three sample types referred to on a given line. For laboratory controls, the code should be CNEG. For any of the test samples provided by SCCWRP, the enter Result.

*Comments*-Place for any explanations that are needed at the observation level.

### ToxWQ

This tab is for the entry of all periodic water quality measurements related to conducting the test.

*SampleID*-Same as for the ToxResults ta.

*ToxBatch*-Same as for the ToxResults tab.

*Dilution*-Same as for the ToxResults tab.

*Concentration*-Same as for the ToxResults tab.

*ConcentrationUnits*-Same as for the ToxResults tab.

*Day*-Day number of the test exposure. For measurements made before test initiation, enter -1. For measurements made on test initiation day, enter 0.

*TimePoint*-Code for which of the two daily measurements the line pertains. For samples that have been in contact with the animals for a day, enter Before. For water that is to be added to the test chambers, enter After. For Days 0 and -1, enter NA.

*Parameter*-Water quality parameter name. Choose from list below.

Code	Water Quality Parameter	Units	Number of Decimal Digits Required
COND	Conductivity	µS/cm	0
DO	Dissolved Oxygen	mg/L	1
PH	pH	pH	2
Hardness	Hardness	mg/L CaCO <sub>3</sub>	0
Alkalinity	Alkalinity	mg/L CaCO <sub>3</sub>	0
AirTemp	Air Temperature	C	1
WaterTemp	Water Temperature	C	1

*Qualifier*-Modifier to the value in the Result field. The expectation for all samples to be tested for this study is that the value should be =, meaning that the value is unmodified. If a situation arises that another modifier is needed, please contact SCCWRP.

*Result*-Value for the analysis for each water quality parameter.

*ResultUnits*-Units for the water quality parameters from the table above.

*LabRep*-Replicate from which water quality sample was taken. If measurement is taken from a surrogate or composite sample, enter 0.

*LabCode*-Same as from the LabTestInfo ta.

*SampleTypeCode*-Same as from the ToxResults tab.

*Comments*- Place for any explanations that are needed for the water quality measurement on the line of data.

### **BroodBoardInfo**

Information on the brood boards used to initiate the tests for the study. Daily neonate counts and health assessment codes are data types of importance must be recorded in this tab.

*BroodBoardID*-Identification number used by the laboratory

*LabCode*-Same as from the LabTestInfo tab.

*DateInitiated*-Date the brood board was started.

*CupNo*-Chamber number for the brood board.

*Day*-Day number of brood board.

*NumNeonates*-Number of neonates for the combination of CupNo and Day.

*HealthCode*-Code for the assessment of the combination of CupNo and Day. Codes from the table below.

Health Code	Description
A	Unhealthy adult
D	Dead adult
K	All OK
M	Male
N	Dead neonates
U	Unhealthy neonates
Y	Neonates used to initiate test

*Comments*-Comments associated with each combination of CupNo and Day.

### **BBtoTestInfo**

This tab has information on brood boards and cup numbers used to initiate the tests.

*BroodBoardID*-Same as from the BroodBoardInfo tab.

*LabCode*-Same as from the LabTestInfo tab.

*ToxBatch*-The test batch identifier for the test that the brood board was used to initiate.

*CupNo*-Cup numbers used to initiate the test.

*Day*-The brood board day when the test was initiated.

### **BroodBoardWQ**

This tab is for water quality data from the brood boards.

*BroodBoardID*-Same as from the BroodBoardInfo tab.

*LabCode*-Same as from the LabTestInfo tab.

*CupNo*- Same as from the BroodBoardInfo tab.

*Day*-Same as from the BroodBoardInfo tab.



*TimePoint*-Same data type as for ToxWQ tab, but for the brood board.

*BBWQParameter*-Same data type as the Parameter from the ToxWQ tab and using the same codes.

*BBWQResult*-Brood board water quality measurement values.

*BBWQUnits*- Same data type as the ResultsUnits from the ToxWQ tab and using the same codes.

*Comments*-Any comments pertaining to the brood board water quality data.

## Appendix D - Study Plan and Logistics

### GENERAL APPROACH

A three-step approach was proposed. During Step 1, all laboratories will participate in an intercomparison exercise by testing a common set of split samples using their current protocols (as described in the current QAPP) and provide more detailed data that may not be routinely collected/reported by all laboratories. Based on the results of Step 1 and discussions between the Expert Science Panel and the Stakeholder Advisory Committee, Step 2 will aim to standardize select *C. dubia* test parameters. Finally, Step 3 will consist of another intercomparison exercise among all laboratories using split samples and the standardized *C. dubia* toxicity testing protocol developed from Step 2.

Appendix D aims to provide an overview of the key study elements and logistics for Step 1 the baseline testing. The specifics of the subsequent steps will depend on the analyses and group discussions of the results of Step 1 baseline testing.

### OVERVIEW OF BASELINE TESTING PROCEDURE

The specific objective of the baseline testing is to collect additional *C. dubia* chronic toxicity data and a more complete/consistent lab technique dataset across California-accredited laboratories. Twelve (12) laboratories will participate in an intercomparison exercise consisting of several split samples tested in three separate testing batches. This testing design is proposed to generate a minimum of seven (7) control datasets per participating laboratory. This was determined statistically based on analyses of the width of the confidence interval to assess intra-laboratory precision. Our analyses indicated that the grand mean for control neonate production from 7 separate tests (each test performed with 10 replicates) would increase our confidence that such mean will fall within the historical control grand mean +/- 5 neonates.

Split samples to be tested include:

Sample 1: Moderately hard dilution water recipe #1 (EPA MH) to be tested at full strength (i.e., 100%). This sample shall be tested along with one (1) laboratory control consisting of the lab's own dilution water recipe.

Sample 2A: Moderately hard dilution water recipe #2 (Perrier) to be tested at full strength (i.e., 100%). This sample shall be tested along with one (1) laboratory control consisting of the lab's own dilution water recipe.

Sample 2B-F: 5 concentrations of sodium chloride (NaCl) diluted in moderately hard water recipe #2 (i.e., Perrier). All samples will be prepared at SCCWRP according to the procedure described earlier in the QAPP. These samples shall be tested as is (i.e., no additional sample dilution allowed) along with one (1) laboratory control consisting of the lab's own dilution water recipe.

Sample 3: NaCl will be provided (as a solid) to each lab with detailed instructions to prepare 5 dilutions using the lab's own dilution water. This serial dilution will be tested along with one (1) laboratory control consisting of the lab's own dilution water recipe. *Note that Sample 3 is now replacing the requirement for each lab to test their routine reference toxicant with each testing batch.*

## SUMMARY OF STANDARD OPERATING PROCEDURES

Participating laboratories (n= 12) will analyze three separate test batches within a ~ 8-week window using their own standard operating procedures for the *C. dubia* chronic toxicity test. A summary of standard operating procedures (SOPs), test acceptability criteria (TAC) and measurement expectations are provided in **Table D1** and in the QAPP. However, all laboratories will be required to meet the following specifications:

- All tests will be carried out to 8 days (i.e., 192 hours).
- All samples, including lab controls, will be performed with 10 replicate chambers.
- Assignment of neonates at test set-up must use the randomized blocking by known parentage, using only brood board chambers with a minimum of 8 neonates from the adult on test initiation day. Each test (i.e., sample and associated laboratory control) will be treated as independent for blocking and randomization except for sample 2A and 2B-F and the two associated controls that must be blocked by the same known parentage.
- A 500 ml-sample of their own dilution water will be collected at test initiation using the container provided by SCCWRP and shipped back to SCCWRP within 24 hours. This sample will be used for analysis of ion composition.
- Test solutions will be renewed daily within a 24 +/- 1 hour window to enhance the comparability of neonate counts among laboratories. Specific time of renewal (hours and minutes) shall be recorded and initialed.

Additionally, participating laboratories will be required to report data that may not be currently documented/reported including (note that the specifics for taking these measurements are provided in the QAPP):

- Number of males, unhealthy and dead adults, and dead neonates in the brood board. This data is to be collected for all days from every chamber within any brood boards that are used to initiate the test. The expectation is that this will be about 6 to 10 days of data depending on the age of the brood board at test initiation
- Specific beginning and end time window for age of neonates at test initiation
- Water quality parameters (air and water temperature, pH, DO, conductivity) at test initiation, termination and before and after daily renewal, to the decimal place specified in the QAPP. If possible, water temperature will also be continuously monitored at the test location.
- Light intensity and twice daily air temperature within the testing area at the time of the experiments and reported in the units specified in the QAPP.

**Table D1.** Summary of test conditions and test acceptability criteria (TAC) for the *Ceriodaphnia dubia* survival and reproduction test.

Parameter <sup>1</sup>	Description
Test organism	<i>Ceriodaphnia dubia</i>
Protocol	EPA/821/R-02-013, EPA 2002 Freshwater Chronic Manual, EPA, 2016 ErrataPA 821-R-02012-ES
Exposure	Static, daily renewal
No. replicate test chambers	10 replicates per sample/dilution
Sample holding time <sup>2</sup>	Up to 48 hours before test initiation
Test duration	8 days, i.e., 192 hours
Endpoints	Survival and reproduction (number of neonates per female)
Laboratory control	One laboratory dilution water control per test sample
Water quality measurements	Daily: air and water temperature in °C, pH and dissolved oxygen in mg/L reported with 0.1 precision; conductivity in µS/cm. Continuous monitoring of water temperature, if possible. Upon receipt and test termination: hardness and alkalinity in mg/L CaCO <sub>3</sub> Once during test in testing area: light intensity in foot-candles; air temperature in °C (0.1 °C precision)
Test Acceptability Criteria (TAC)	80% or greater survival and an average of 15 or more live neonates per surviving female in the controls at test termination (i.e., 8 days)

<sup>1</sup> Parameters and test conditions used in this study are suitable for investigative/non-compliance testing but may be different than those required for NPDES permit testing.

<sup>2</sup> This is a deviation from the promulgated method.

## OVERVIEW OF SPLIT SAMPLE PREPARATION AND DISTRIBUTION

Bulk test water samples will be prepared in the SCCWRP laboratories as described in the QAPP using large sample containers with spigots and thoroughly mixed on a large-capacity stirrer to ensure that the samples are homogeneous. The number of samples to be tested by the participating laboratories are presented in **Table D2**. Bulk samples will be allowed to equilibrate for up to 48 hours before preparing the split-samples that will be shipped to the laboratories. Subsampling of the bulk test samples will be conducted using 3.8 L cubitainers filled to the top. All cubitainers will be filled at random in two steps. First, each cubitainer will be filled halfway. Then the cubitainers will be filled the rest of the way in no particular order. Each cubitainer will be labelled with a unique sample ID and stored in the dark in the walk-in fridge at 4 °C less than 48 hours before shipping them to the participating laboratories.

**Table D2.** Number of split-samples to be tested by the 12 participating laboratories for each round. Three testing rounds will be completed for this study.

Sample ID	Number of samples per lab per round	Number of sample dilutions to test	Number of lab control to include per sample
1	1	_*	1
2A	1	_*	1
2B-G	5	_*	1
3	1 <sup>‡</sup>	5	1

\*Water samples DO NOT require further dilution before testing.

<sup>‡</sup> Sample 3 will be shipped as a powder with instructions to prepare the serial dilution for testing.

To evaluate their preparation method and prevent unexpected toxicity, SCCWRP will prepare bulk water samples and send them to one laboratory for a *C. dubia* chronic toxicity test. If unspiked samples are not toxic and a dose-response appears normal for the dilution series, the preparation method will be deemed suitable for the ILS. SCCWRP will then prepare fresh bulk samples and split them in individual cubitainers as described above. Since all methods and equipment will be the same for subsequent rounds, this preliminary testing will only be carried out for round one.

To ensure that all subsamples are representative of the original bulk test samples, two subsamples will be collected in separate vessels from each cubitainer before shipment. The first set (50 ml) will be used to measure conductivity, alkalinity, and hardness. The subsample must be discarded after the measurements are completed. The second set (500 mL) will be collected for ion composition analysis. Due to sample volume requirements, ion composition analysis will only be collected by SCCWRP from each cubitainer before shipping. These subsamples will be collected in 500 ml HDPE bottles, filled to the top, and shipped to the analytical laboratory (Physis) to measure bicarbonate, carbonate, chloride, fluoride, nitrate, sulfate, selenium, and major cations (calcium, phosphate, magnesium, sodium). The analyses will be completed within 14 days of sampling to meet holding time requirements.

Split samples will be shipped to each laboratory starting August 22, 2022 according to the schedule presented below. Samples will be shipped on wet ice using priority overnight (OnTrac or FedEx) service to the laboratories to the addresses in Appendix A of the QAPP. The shipments will also include chain-of-custody (COC) forms completed by SCCWRP and a copy of the study plan and testing instructions. SCCWRP will notify the laboratories via email once the samples are in transit and provide a tracking number. It is the responsibility of the laboratories to contact SCCWRP if they have not received the samples by the following day 2:00 pm.

Upon delivery, temperature, conductivity, pH, dissolved oxygen, hardness, and alkalinity, must be measured and recorded for each sample to document their stability before testing is initiated. These measurements shall be made from a small subsample poured into a clean secondary vessel. Probes and any other measuring equipment cannot be used in the cubitainer, and the subsample used for water quality must be discarded after use (subsample cannot be used for testing or as a chemistry or archived sample). Additionally, a 125-ml sample must be collected from each cubitainer at the time of test initiation and archived. Once all chemistry and water quality samples have been collected by both SCCWRP and the laboratories, there should be more than 3 L remaining in each cubitainer to conduct the 8-day *C. dubia* test.

For sample 3, 7.00 g of NaCl will be weighed and placed in 100 ml HDPE containers. Each laboratory will receive one container with instructions to prepare the serial dilution using their own lab dilution water (i.e.,

dilute the supplied NaCl in 3.5 L of their own dilution water and perform a 50% dilution series to generate a total of 5 dilutions). Similar to the split-water samples, once the dilutions are prepared, the laboratories must record temperature, conductivity, pH, dissolved oxygen, hardness, and alkalinity for each dilution at test initiation. A 125-ml sample must also be collected and archived from each dilution at test initiation.

Note that approximately one (1) hour prior to test initiation and water changes, the volume of water needed to renew the test solutions should be adjusted/maintained at test temperature.

## **DATA SUBMISSION**

SCCWRP will provide an Excel data submittal form and culture/bench sheet templates to the participating laboratories. All test data in electronic format and scanned copies of the culture/bench sheets must be submitted to the SCCWRP data portal no later than October 18. Data required include:

- Laboratory information
- Sample information upon receipt (time, temperature, condition, and more as described above)
- Testing conditions including dilution water and food recipe
- Brood board health data
- Bench water quality data for testing, survival and reproduction counts
- Control charts for reference toxicant tests for the last 12 months

Detailed data requirements are provided in the QAPP. Note that continuous water temperature data shall be submitted separately by email to SCCWRP.

## **COMMUNICATION**

Participating laboratories and other stakeholders will meet with SCCWRP, and the Expert Science Panel advising on this project to finalize the study plan, discuss logistics and review the results. A minimum of three remote meetings will be scheduled to provide a forum for discussion and clear communication among the project team and participants. Additional communication via email will be encouraged throughout the study. For more information on the overall study design and coordination meetings, please contact Alvina Mehinto [alvinam@sccwrp.org](mailto:alvinam@sccwrp.org). For questions regarding samples shipping from and to SCCWRP and data submission, please contact Darrin Greenstein [darring@sccwrp.org](mailto:darring@sccwrp.org).

The first meeting, held remotely on May 24, 2022, and attended by the stakeholders aimed to review the first draft of the testing approach (including sample preparation and shipping, test measurements and data reporting) and discuss the timeline for testing and data submission. The second meeting held on June 24, 2022, among members of the Expert Science Panel, stakeholders and laboratories aimed to further refine the study design. A third meeting, held on July 11, 2022, as a closed session at the request of the Expert Science Panel, aimed to finalize the testing design. The fourth meeting, held on August 3, 2022, aimed to train the participating in data collection and data submission, and answer logistics questions.

## SCHEDULE

- ~~May 17: Draft study plan sent to all stakeholders for review~~
- ~~May 24: Stakeholder Committee meeting, held via Zoom, to discuss the first draft of the study plan~~
- ~~June 14: Revised study plan submitted to the Expert Science Panel~~
- ~~June 24: Public meeting with Expert Science Panel and stakeholders' representatives to refine the study plan~~
- ~~June 11: Expert Panel closed session to finalize the study design~~
- ~~July 19: Revised study plan and draft QAPP sent to stakeholders and Science Panel. Final comments are due July 28 at 5pm PDT.~~
- ~~August 2: Meeting with participating laboratories. SCCWRP review testing requirements and provided training for data collection/submission.~~
- ~~August 3: Revised QAPP submitted to the State Water Board~~
- ~~August 8: Approval of the QAPP by the Expert Science Panel and the State Waterboard~~
- **August 12-13:** A batch of split samples prepared by SCCWRP for preliminary test
- **August 15:** Samples sent to selected laboratory for preliminary testing
- **August 18-19:** First batch of split sample prepared by SCCWRP
- **August 22:** Cubitainers containing first batch of split samples shipped to the laboratories.
- **August 23:** First batch of *C. dubia* toxicity tests performed
- **September 9-10:** Second batch of split samples prepared by SCCWRP
- **September 12:** Cubitainers containing second batch of split samples shipped to the laboratories
- **September 13:** Second batch of *C. dubia* toxicity tests performed
- **September 23-24:** Third batch of split samples prepared by SCCWRP
- **September 26:** Cubitainers containing third batch of split samples shipped to the laboratories.
- **September 27:** Third batch of *C. dubia* toxicity tests performed
- **October 18:** Deadline for data submission

## CONTINGENCIES

**Lost Samples:** If a sample is not delivered to a laboratory on the expected arrival date or if the sample has spilled during shipment, the laboratory must contact SCCWRP promptly. SCCWRP will ship new cubitainers that same day. However, this second batch of samples sent must be tested within 24 hours to ensure that holding times are comparable to other laboratories.

**Failed Test Acceptability Criteria (TAC):** A laboratory will be able to retest up to two test batches if the test conditions and TAC are not met (see **Table 1** and the QAPP for this study). A laboratory planning to retest must contact SCCWRP within 24 hours of knowing that a test failed the TAC. Laboratories are encouraged to retest with the remaining sample; however, arrangements might be made to re-test with archived samples. Laboratories that fail to provide data for Step 1 baseline testing may still be considered to participate in the confirmation testing in Step 3.

**Late Data Submission:** All data must be submitted to the SCCWRP data portal and pass the QA checkers by **October 18, 2022**. If a laboratory experiences some delays, SCCWRP must be contacted no later than 48 hours before the deadline. Laboratories will be granted an additional three (3) days to submit all their data. Past this new deadline, SCCWRP cannot guarantee that the data will be used in subsequent data analyses.