EPA METHOD 314.0 PERCHLORATE

LA	B NAME:DATE:			
1.	1. Columns and Reagents Sections 3, 6, 7			
Γ	Holding Time 28 days			
Γ	Column AS 16			
Γ	Volume of sample 1.0 ml			
Γ	Eluent – 50 mM sodium hydroxide			
Γ	Matrix Pretreatment cartridges – Ba (SO ₄), Ag (CI), H (cations and Ag)			
Γ	Mixed Common Anion Stock – Cl/SO₄/CO₃ 25mg/ml each			
Γ	KCl – 1410 μ S/cm at 25°C ± 30 μ S/cm			
Γ	Std Perchlorate expires 1 yr COMMENTS:			
2.	Quality Control			
Γ	Initial Demonstration Section 9			
Γ	Low Background			
Γ	7 replicates of 25.0 μ g/L. % R = 80 –120% of TV % RSD = < 10%			
Γ	Second Source			
Γ	MDL – initially determined when a new operator begins work, significant change in background or instrument response. DO NOT SUBTRACT BLANK VALUES!!			
Γ	MRL – Minimum Reporting Level - > 3X MDL or concentration, which would yield a S/N of 5. MRL must never be lower than the lowest calibration standard			
Γ	MCT – Matrix Conductivity Threshold			
Γ	If MRL is between 3 – 5 ppb use 25 ppb Perchlorate.			
Γ	Prepare a series of solutions with constant Perchlorate conc. And increasing common anions (200, 300, 400, 500 1000 mg/L).			
Γ	Reference 400 mg/L conductance between 3200 –3700 uS/cm. LFB/200 mg/L mixed common ion should have same A/H and CIO_4 conc.			
Γ	Determine % difference (PD) between the average A/H ratios for LFB and each mixed common anion solutions.			
Γ	MCT = the matrix conductance where the PD exceeds 20%. Plot PD vs. matrix conductance. $r^2 > 0.95$			
	Γ MCT = 20% X (X coefficient) +(constant intercept value)			
	Γ Alternately, MCT conductance level of the highest mixed anion solution which yielded a PD value below the 20% threshold.			
	Γ MRL confirmation % R = 70 – 130%			

	Spike soln. at MRL with common anions a conc. est. to generate a conductance near the MCT within \pm 10%).		
_	If fails criteria, lower the MCT by 10% and MRL re-verified.		
]	LRB = < ½ MRL		
]	IPC (same soln. as MRL confirmation) initially PD = < 25% Repeat MCT if any one fails		
]	LFB (use same std as calibration) %Rec. = 85 – 115%		
]	LFM (spike conc. = or > spl. conc.) %Rec. = 80 - 120%		
]	Dup or LFM dup RPD = $\pm 15\%$		
(COMMENTS:		
3	3. Calibration Section 10		
Γ	Initial 5 points for concentrations from 4 –400 ppb. Linearity assumed if % RSD < 15% average ration or RF can be used. Use Peak Areas.		
Γ	Continuing ICC – use lowest standard %Rec. = 75 – 125% CCS after every 10 and end %Rec. = 85 – 115%		
(COMMENTS:		
4. Procedure Section 11			
Γ	Conductance of all samples measured. If above the MCT – dilution or pretreatment is performed.		
Γ	Dilution done – Check if MRL raised by a proportion equivalent to the dilution.		
Γ	If Perchlorate found, back calculate. If none found and monitoring is required below the MRL, pretreatment is required.		
Γ	Pretreatment using BA/Ag/H cartridges same QC limits applied.		
Γ	LRB and LFB and LFM must be pretreated with any pretreated samples.		
(COMMENTS:		
	5. Data Analysis and Calculations Section 12		
Γ	RT window within ± 5%		
Γ	Concentration based on ICAL		
Γ	Report values between MRL and highest standard.		
Γ	Reported in ug/L		

COMMENTS:	
6. Notes Section 1, 9, 11	
Γ For any changes, Section 9.2; 1.2.1 and 11.2.6 must be adhered to. Γ Maintenance is required when there is reduced peak response and peaks are non-symmetrical.	
COMMENTS:	