COVID-19 Surveillance Update

Presentation to SCCWRP Commission

September 3, 2021

Joshua A. Steele

Background

- COVID-19 pandemic has put wastewater based surveillance of viral pathogens at the forefront of public health
- California's wastewater community continues to be active in wastewater based surveillance for SARS-CoV-2
 - More than 40 facilities are monitoring and have at least six months of data
 - Several facilities (including SCCWRP member POTWs) have more than a year of data
 - Participation in national studies
- WQMC Wastewater Based Epidemiology continuing to provide guidance on effective data use and interpretation


From Method Comparison to Effective Data Use

Laboratory measurement method variability

- Measurement methods sensitivity
 - Guidance document with WBE Subcommittee

SARS-CoV-2 variants in wastewater

Potential Sources of Variability



Manuscript in revision

From Method Comparison to Effective Data Use

Laboratory measurement method variability

- Measurement methods sensitivity
 - Guidance document with WBE Subcommittee

SARS-CoV-2 variants in wastewater

Recommendations to Improve Sensitivity


- Remove heat inactivation step
 - 50% of treatment plants using heat inactivation
- Concentrate or measure a larger volume
 - ~30% of labs are not concentrating
- Measure SARS-CoV-2 using digital RT-PCR
 - 60% of labs are using RT-QPCR
- Use sludge samples
 - 25% of labs are collecting sludge


Sensitivity Case Study at JWPCP

Opportunity to test sensitivity recommendations

- Split samples among four groups
 - Small volume qPCR (Zymo Research)
 - Small volume, concentrated qPCR (BioBot)
 - Large volume qPCR (UC Berkeley COVID-WEB)
 - Large volume digital PCR (SCCWRP)

Preliminary results





- UC Berkeley and SCCWRP are highly correlated (r²=0.7)
- 60% of BioBot samples below detection

 50% of Zymo samples are below detection

 Zymo's concentrations are much higher compared to others

Takeaways from the Sensitivity Comparison

 Larger volume concentration qPCR and digital PCR are sensitive enough at low concentrations

 Smaller volume qPCR methods too insensitive to be useful when prevalence is low

WBE Subcommittee recommendations looking pretty good

From Method Comparison to Effective Data Use

Laboratory measurement method variability

- Measurement methods sensitivity
 - Guidance document with WBE Subcommittee

SARS-CoV-2 variants in wastewater

Focus on SARS-CoV-2 Variant Analysis

- Variant analysis is the next on the horizon for SARS-COV-2 analysis
 - High public health impact
 - Not a part of current SARS-CoV-2 monitoring
- Should variant analysis become a part of routine wastewater based surveillance?
- WBE Subcommittee has determined that variant analysis is not ready to be included in routine monitoring
 - Guidance document sent to CASA and Public health departments.

SARS-CoV-2 Variant Analysis in Sewage

- Wastewater is complex and multiple variants in sewage
 - Detection of a single mutation is not enough
- Approaches using both sequencing and targeted analysis are needed
 - Active area of research
- Interpretation of variant data in sewage will require collaboration across multiple agencies

Ongoing Variant Research at SCCWRP

- Sequencing SARS-CoV-2 genomes in sewage samples
 - Sewage RNA Virus manuscript in revision
 - Identifying SARS-CoV-2 mutations in sewage
- Quantifying Variants using digital PCR
 - Targeting prevalent variants
 - Using multiple gene targets to identify specific variants

 Collaborating with universities, POTWs, and public health agencies within California and across the US